

Product datasheet for **MR216079**

Ckap2 (NM_001004140) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ckap2 (NM_001004140) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ckap2
Synonyms:	AI461656; AW228814; LB1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR216079 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC

ATGGCAGAGTCCAGGAAACGCTTCTTGGCAGGGCTGCGAGGAATCCCTTGCCGGTGACCCGGGACCTGC
AGCTGCCCCCGACCAGGCGAGACCAGCCGGCCTTCCGAGAACAAGAAAAACAAAACTCAAGGAATATCT
GTTAATACGAAAAACAGTGTTTCCATACAAACAAGAAAATCAGATATCCAGAGATCAGAAAATGATAACA
TCTGAGGATCGGGTCCAAGAAGGAAAAAAGTTGTGAAACTCAAAACAGAAGTGGCTGATAAAGAAAAATA
TTGAAAGTACAGTTGAGAAAACTGCATTCCATTAAGCTGGAGAAGTAACCAGCTCAGAAAATACATAA
TTCAAAGGATAATGTTCAAGCTGTGCAGTTATTATCTACAAGAGATGATTTGCCAGGTGACTGTGACA
TTAGACCCGGCATGTACCATAAAGACAATAAAAAAGATGCAAATGACTGCAGAGAAGCCCAACAAGACA
GTAATGTGTCTAAAAAGCGAGTGTGGCTATTATCATGGCCAGATTGTTCACTAAGATTAATTCATT
TAGAAAGCTTCCAGTGTCAAAGGCGAGAGTCTTACAACAACCAAGAACTTCCCACTACTGTGTCTAAA
GCCATGAAGGCTCAGAGTGAGCCTGCGAACACCGTCAGCGTAAAAGCCTCAACTACAGCTGCTGCCACCA
AGTTTGCAGACGCTAAACCTGTGAGCACTGCGTCGAAGGACACACTTGTGCGACCACCTATCCGAAGCCT
CCACAGTAGTTCCACGGAGCTGCGAAACAGGGCCTCAGCAGGCCGCTTGCCAATGTTACAGTCAGAAAA
GGGATGCTTGACAAGGAGTCACATCGATCAGAACCAGTTGTGTCGGTGTCAAAGCCGGTCTTCTCAGG
CCCCATCAAGAAGCATAGCATCTAAGATGCAGCCAGGACTGATTCGTCTAATACCAGACTGATGGTAAA
GCCAAAAGATACTGACCAGCGCAGATACACGATTGCAGGAGCAAGTTCATAGATCAGCTCAGCTCAAG
GACACCACCGCAGAGAGGAAAGCTCAGATGACGGAGTGGAGAAGTGGCAAAGGAAAGGGGCTAAAACGAC
CTCCTCATTCAGTAGTCACCCAGGCTGAGCCAAAGGGACAAAAGTAAAACCCAGTTGGGTCCTTTTGGAC
TACCATTGGCAGAAGAAGATGAGCAGAGATTATTTACTGAAAAAGTAAAACAAAACTATTTCTGAATGCCTG
AACCTTATTAATGAGGGATGCCCAAAAGAAGAAATTTTGGCTACATTGAATGACCTGATTCATAATATTC
CAGATGCCAAAAAAGTGTAAATATTGGATATGTCTTGTACGATTGAACCAATCACAAGTCTATTGA
AAATATTATCTCAATCTATGAGAAGGCCATCCTGGCAGGGGCTCAGCCTATTGAAGAGATGCGTCATATA
ATCATAGATATTCTAACACGAAGAGTCAAGAAAAAGTCAATTTGGGAGAAAATATTGAGGAGGCTCATG
CAACCAAGGAACCTATCCAAGAAGTAAATGCTGATGCTAATGTAGGGTCAGGAAAACAGGAGAGGAGAA
TGAACATCATGGCAAAGTTGAGGTGTATGAAGATGATCAAGACAACAAAAATAAAGATCCGGATTTAACA
ACTCCGGACTCAAAAAGTGAAGCAGGTTGCATAATTCGATATAATGTATCTTCTACACCTCGGCTGCAAA
GTATGAAAAAGATGCAACATGATAAAAAATCCACACTTAAAGAGCTAAAGTTCTGACACCAAGTGAGACG
CTCACGGCGCATTCAAGACAAGACTTACGACTGCCAGCCATGTTAAAAGACCACGACCTTCTGTGTCT
TCCTGGAGCAGTTGTCTGAGTTGGGAGGGGATGCCTTTGTCTGCCGCCCTAATGCAGCACTGTGCCCTC
TGTTCTTTGAGACCGACGTAGCGGAAGAGGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR216079 protein sequence
Red=Cloning site Green=Tags(s)

```


MAESRKRF LGRAARNPLPVTRDLQLPPTRRDQPAFREQRKQKLKEYLLIRKTVFPYKQENQISR DQKMIT
SEDRVQEGKVVV LKTEVADKENIESTVEKNCIPLKAGEVTSSEIHNSKDNVQAVQLLSTRDDLPGQTVT
LDPACHHKDNKMQMTAEKPKQDSNVSKRVLGYHGGQIVQSKINSFRKLP SVKGESLTTTKLPTTVSK
AMKAQSEPANTVSVKASTTAAATKFADAKPVSTASKDTLVRPPIRSLHSSSHGAAKQGLSRPLANVTVRK
GMLDKESHRS EPVSVVKAGSSQAPSRSIASKDAARTDSSNTRL MVKPKD TDQRRYTIAGASVHRS AQLK
DTTAERKAQMTEWRTGKGKGLKRPPHSVVTQAEPKGQSEN PVGSFWTTMAEEDQRLFTEKVNKTI SECL
NLINEGCPKEEILATLNDLIHNIPDAKKLVKYWICLVRIEPI TSP IENIISIYEKAILAGAP IEE MRHI
IIDILTTKSQEKVNLGENIEEAHATKEPIQEVNADANVGS GKPGEENEHGHKVEVYEDDQDNKIKDPDLT
TPDSKTEAGCIIRYNV SSTPRLQSMKMQHDKNSTL KELKFLTPVRRSRRIQDKTSRLPAMLKDHDPSVS
SLEQLSELGGDAFVCRPNAALCPLFFETDVAEEE
    
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001004140

ORF Size: 1995 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001004140.2](#), [NP_001004140.2](#)

RefSeq Size: 2601 bp

RefSeq ORF: 1995 bp

Locus ID: 80986


UniProt ID: [Q3V1H1](#)

Cytogenetics: 8 A2

MW: 74 kDa

Gene Summary: Possesses microtubule stabilizing properties. Involved in regulating aneuploidy, cell cycling, and cell death in a p53-dependent manner.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR216079