

## Product datasheet for **MR216103**

### **LnX1 (NM\_001159579) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	LnX1 (NM_001159579) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	LnX1
Synonyms:	LnX
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**ORF Nucleotide Sequence:**

>MR216103 representing NM\_001159579  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGAAGGCGCTGCTGCTTCTGGTCTGCCTGGCTGAGTCCAGCTAACTACATTGACAATGTGGCAACC  
 TGCACTTCCTGTACTCGAACTCTGCAAGGGGCTCCCACTATGGCCTGACCAAGGACAGGAAAAGCGG  
 CTCTCAAGATGGCTGCTCTGATGGTTGCGCCAGCCTCATGGCAACGACGCTCTCCCCAGAGGTTTCTGCA  
 GCCGCCACCATCTCTTTAATGACAGATGAGCCTGGCCTAGACAACCCTGCCTACGTGTCCTCGTGGAGG  
 ATGGTGAGCCAGTAGCCAACCTCAGACTCTGGCCGGAGCAACCGGACTAGGGCACGGCCCTTTGAGCG  
 CTCCACTATGAGAAGCCGGTCTTTCAAGAAGATCAACCGAGCTCTAAGTCTCTCCGGAGGACGAAGAGC  
 GGGAGTGTGTTGCCAACATGTTGACCAGGGCAGGGACAACCTGAGAACACCCTGTCCCGAAGTCT  
 TCCAAGTGTGTTTACCTGATCCAGATGGTAAATACAAGCATCAAATCAACAGAGCGGATCCAG  
 CGAAAGCCTCTCATAAGGCTGGTGGGGGCGAGGAGACGCCGCTGGTCCACATCATCCAGCATATT  
 TATCGCATGGAGTGATTGCCAGAGATGGACGGCTGCTGCCAGGAGACATCATTCTCAAGGTCAACGGCA  
 TGGACATCAGCAATGTCCCTCATAACTATGCCGTGCGGCTCCTGCGGCAGCCCTGCCAGGTGCTGCGGCT  
 AACGGTGTGCGGGAGCAGAAGTTCCGTAGCAGAAGCAATGCGCACGTGCCGGACTCCTATGGACCTCGG  
 GATGACAGCTTCATGTAATCTCAACAAAAGCAGCCCGAGGAGCAGCTTGGCATAAAGCTGGTGGCC  
 GGGTGGATGAGCCCGGCTGTTTCATCTCAACGTGCTAAATGGGGTGTGCGCGACCGACACGGCCAGCT  
 GGAGGAGAATGACCGTGTGCTGGCCATCAACGGCATGACCTTCGATTTGGCAGTCCAGAGAGTGACGCT  
 CATCTGATTCAGGCCAGTGAAAGACGCGTCCATCTCGTCGTGTCCCGCAGGTTGACAGTCCAGCCAG  
 ACATCTTTCAAGAAGCTGGCTGGATCAGCAATGGCCAGCAATCCCAGGCCAGGAGAGAGGAACACAGC  
 TTCTAAGCCTGCAGCCACCTGTGATGAGAAGTTGTAAGTGTCTGGAAAGACCCAGCGAGTCTCTCGGC  
 ATGACCGTGGAGGGGAGCATCCACAGGGAATGGGACTTGCCATCTATGTCATCAGTGTGAACCCG  
 GCGGAGTCATAAGCAGAGACGGGAGAATAAAAAACAGGTGACATTTTATTGAATGTGAATGGGATTGAGCT  
 CACAGAGGTCAGCCGGACAGAGGCCGTTGCCATACTGAAGAGCGCACCCCTCCTCAGTGGTACTCAAAGCC  
 TTGGAAGTCAAGGAACAGGAGGCCAGGAAGACTGCAGCCAGCAGCCCTGGACTCCAACCACAACGTGA  
 CCCCCTGGTACTGGTCCCCATCTGGGTGATGTGGCTGGAATTACCACAGTACTTATGTAAGTGA  
 AGATGTGATACTGCGAAGGAACACAGCTGGAAGCCTGGGCTTCTGCATTGTGGGAGTTATGAAGATAC  
 AGTGGGAACAACCTTTTTTATCAAGTCCATTGTTGAAGGAACCTGCATACAATGACGGAAGATCA  
 GATGTGGTACATTCTTCTCGCTGTCAACGGTAGAAGTACATCGGGTATGATACAGCTTGCCTGCCAG  
 GATGCTCAAGGAACCTAAAGGGAGAATTACTCTGACCATTGCTTCTTGGCCTGGTACTTTTTTA

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
 TGGATTACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:**


>MR216103 representing NM\_001159579  
 Red=Cloning site Green=Tags(s)

MKALLLLVLPWLSANYIDNVGNLHFLYSELCKGASHYGLTKDRKRSQDGPDPGCASLMATTLSPVEVA  
 AATISLMTDEPGLDNPAYVSSVEDGEPVANSSDSGRSNRTRARPFERSTMRSRSFKKINRAL SALRRTKS  
 GSVVANHVDQGRDNSENTTVPEVFPRLFHLIPDGEITSIKINRADPSELSIRLVGGSETPLVHIIQHI  
 YRDGVIARDGRLLPGDIILKVNMGMDISNVPHNYAVRLLRQPCQVLRRLTVLREQKFRSRSNAHVPDSYGPR  
 DDSFHVILNKSSPEEQLGIKLVRRVDEPGVIFNVLNNGVADRHGQLEENDRVLAINGHDLRFGSPESAA  
 HLIQASERRVHLVVSQRVQSSPDIFQEAGWISNGQQSPGGERNTASKPAATCHEKVVSVMKDPSESLG  
 MTVGGGASHREWDLPIYVISVEPGGVISRDGRIKTDGILLNVNGIELTEVSRTEAVAILKSAPSSVVLKA  
 LEVKEQEAQEDCSPAALDSNHNVTPPGDWSPSWMWLELPQYLCNCKDVILRRNTAGSLGFCIVGGYEEY  
 SGNKPFPIKSIVEGTPAYNDGRIRCGDILLAVNGRSTSGMIHAACLARMLKELKGRITLTIASWPGTFL


SGP**TRTRRLEQKLI**SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001159579

ORF Size: 1884 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001159579.1</a> , <a href="#">NP_001153051.1</a>
<b>RefSeq Size:</b>	2347 bp
<b>RefSeq ORF:</b>	1887 bp
<b>Locus ID:</b>	16924
<b>UniProt ID:</b>	<a href="#">O70263</a>
<b>Cytogenetics:</b>	5 C3.3
<b>MW:</b>	69.2 kDa
<b>Gene Summary:</b>	E3 ubiquitin-protein ligase that mediates ubiquitination and subsequent proteasomal degradation of NUMB. E3 ubiquitin ligases accept ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates. Mediates ubiquitination of isoform p66 and isoform p72 of NUMB, but not that of isoform p71 or isoform p65.[UniProtKB/Swiss-Prot Function]