

Product datasheet for **MR216141L3V**

Enpep (NM_007934) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Enpep (NM_007934) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Enpep
Synonyms:	6030431M22Rik; APA; Bp-1/6C3; Ly-51; Ly51
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_007934
ORF Size:	2835 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216141).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007934.3 , NP_031960.1
RefSeq Size:	4208 bp
RefSeq ORF:	2838 bp
Locus ID:	13809
UniProt ID:	P16406
Cytogenetics:	3 57.92 cM
Gene Summary:	Regulates central hypertension through its calcium-modulated preference to cleave N-terminal acidic residues from peptides such as angiotensin II.[UniProtKB/Swiss-Prot Function]



[View online »](#)