

Product datasheet for **MR216146L3V**

Acbd5 (NM_001102437) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Acbd5 (NM_001102437) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Acbd5
Synonyms:	1300014E15Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001102437
ORF Size:	1527 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216146).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001102437.1 , NP_001095907.1
RefSeq Size:	3746 bp
RefSeq ORF:	1530 bp
Locus ID:	74159
UniProt ID:	Q5XG73
Cytogenetics:	2 A3
Gene Summary:	Acyl-CoA binding protein which acts as the peroxisome receptor for pexophagy but is dispensable for aggrephagy and nonselective autophagy. Binds medium- and long-chain acyl-CoA esters (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)