

## Product datasheet for **MR216173L3V**

### Ap4b1 (NM\_026193) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ap4b1 (NM_026193) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ap4b1
Synonyms:	1810038H16Rik; Ap4b4; AV004952
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_026193
ORF Size:	2217 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216173).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_026193.2</a> , <a href="#">NP_080469.2</a>
RefSeq Size:	2914 bp
RefSeq ORF:	2217 bp
Locus ID:	67489
UniProt ID:	<a href="#">Q9WV76</a>
Cytogenetics:	3 F2.2



[View online »](#)

**Gene Summary:**

Component of the adaptor protein complex 4 (AP-4). Adaptor protein complexes are vesicle coat components involved both in vesicle formation and cargo selection. They control the vesicular transport of proteins in different trafficking pathways. AP-4 forms a non clathrin-associated coat on vesicles departing the trans-Golgi network (TGN) and may be involved in the targeting of proteins from the trans-Golgi network (TGN) to the endosomal-lysosomal system (By similarity). It is also involved in protein sorting to the basolateral membrane in epithelial cells and the proper asymmetric localization of somatodendritic proteins in neurons (PubMed:18341993). AP-4 is involved in the recognition and binding of tyrosine-based sorting signals found in the cytoplasmic part of cargos, but may also recognize other types of sorting signal (By similarity).[UniProtKB/Swiss-Prot Function]