

Product datasheet for **MR216176L4V**

Tymp (NM_138302) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tymp (NM_138302) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tymp
Synonyms:	2900072D10Rik; Ecgf1; PD-ECGF; PDECGF; Pdgfec
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_138302
ORF Size:	1413 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216176).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_138302.1
RefSeq Size:	2042 bp
RefSeq ORF:	1416 bp
Locus ID:	72962
UniProt ID:	Q99N42
Cytogenetics:	15 44.86 cM
Gene Summary:	Catalyzes the reversible phosphorolysis of thymidine. The produced molecules are then utilized as carbon and energy sources or in the rescue of pyrimidine bases for nucleotide synthesis (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)