

Product datasheet for MR216179

Impmp1l (NM_028260) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Impmp1l (NM_028260) Mouse Tagged ORF Clone
 Tag: Myc-DDK
 Symbol: Impmp1l
 Synonyms: 1500034J20Rik; 2610528O17Rik; AU042386
 Mammalian Cell Selection: Neomycin
 Vector: pCMV6-Entry (PS100001)
 E. coli Selection: Kanamycin (25 ug/mL)
 ORF Nucleotide Sequence: >MR216179 representing NM_028260
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGGATCGCC

ATGCTTCGTGGTGTCTGGGAAAGGCTTTTCGACTTGCTGGCTATACCATTCAGTATGGCTGCATAGCTC
 ATTGTGCTTTTGAATACGTTGGTGGTGTGTAATGTGTTCTGGACCATCAATGGAACCCACAATCAAAA
 TTCAGATATTGTCTTTGCAGAAAATCTTAGCCGACATTTTATGGCATCCAAAGAGGTGACATTGTAATT
 GCAAAAAGCCCAAGTGATCCAAAATCAAATATCTGTAAGAGTAATTGGTTTGAAGGAGACAAAATCC
 TCAGCACTAGTCCATCTGATGTATTTAAGAGCCGAGTTACGTGCCAACAGGTCATGTTTGGTTAGAAGG
 TGATAACTTACAGAATTCTACCGATTCCAGGTATTATGGACCTATACCATATGGACTAATAAGAGGACGT
 ATCTTCTTTAAGATTTGGCCTTTCAGTGACTTCGGATTTCTGCGTGACAGCCCAAATGGCCAGAGATTTT
 CTGATGAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/mm9051_d02.zip

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_028260

ORF Size: 498 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_028260.2](#), [NP_082536.1](#)

RefSeq Size: 955 bp

RefSeq ORF: 501 bp

Locus ID: 66541


UniProt ID: [Q9CQU8](#)

Cytogenetics: 2 E3

MW: 19 kDa

Gene Summary: Catalyzes the removal of transit peptides required for the targeting of proteins from the mitochondrial matrix, across the inner membrane, into the inter-membrane space. Known to process the nuclear encoded protein DIABLO (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR216179