

Product datasheet for **MR216203L3V**

Mesdc2 (Mesd) (NM_023403) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mesdc2 (Mesd) (NM_023403) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mesd
Synonyms:	2210015O11Rik; AW537813; mesd; mKIAA0081; msd
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_023403
ORF Size:	672 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216203).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_023403.3
RefSeq Size:	4492 bp
RefSeq ORF:	675 bp
Locus ID:	67943
UniProt ID:	Q9ERE7
Cytogenetics:	7 D3



[View online »](#)

Gene Summary:

Chaperone specifically assisting the folding of beta-propeller/EGF modules within the family of low-density lipoprotein receptors (LDLRs). Acts as a modulator of the Wnt pathway through chaperoning the coreceptors of the canonical Wnt pathway, LRP5 and LRP6, to the plasma membrane. Essential for specification of embryonic polarity and mesoderm induction (PubMed:12581525). Plays an essential role in neuromuscular junction (NMJ) formation by promoting cell-surface expression of LRP4 (PubMed:24140340). May regulate phagocytosis of apoptotic retinal pigment epithelium (RPE) cells (PubMed:27184668).[UniProtKB/Swiss-Prot Function]