

Product datasheet for **MR216205L4V**

Fbxl20 (NM_028149) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Fbxl20
Synonyms:	2610511F20Rik; 4632423N09Rik; A1849362; A117906; C86145; Fbl2; mKIAA4147; Scr
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_028149
ORF Size:	1308 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR216205).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_028149.1](#), [NP_082425.1](#)

RefSeq Size: 8607 bp

RefSeq ORF: 1311 bp

Locus ID: 72194

UniProt ID: [Q9CZV8](#)

Cytogenetics: 11 D



Gene Summary:

Substrate-recognition component of the SCF (SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex. Isoform 3 regulates neural transmission by binding and ubiquitinating RIMS1, a modulator of presynaptic plasticity.[UniProtKB/Swiss-Prot Function]