

Product datasheet for MR216238L4

Arhgef6 (NM_152801) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Arhgef6 (NM_152801) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Arhgef6
Synonyms:	1600028C08Rik; 1700038J06Rik; 4930592P22Rik; alpha-PIX
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216238).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	NM_152801
ORF Size:	2385 bp



[View online »](#)

