

Product datasheet for **MR216258L4V**

Etv2 (NM_007959) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Etv2 (NM_007959) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Etv2
Synonyms:	Etsrp71
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_007959
ORF Size:	1005 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216258).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007959.2 , NP_031985.2
RefSeq Size:	1050 bp
RefSeq ORF:	1008 bp
Locus ID:	14008
UniProt ID:	P41163
Cytogenetics:	7 B1
Gene Summary:	Binds to DNA sequences containing the consensus pentanucleotide 5'-CGGA[AT]-3'. [UniProtKB/Swiss-Prot Function]



[View online »](#)