

## Product datasheet for **MR216265L3V**

### **Kmt5b (NM\_001167885) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Kmt5b (NM_001167885) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Kmt5b
Synonyms:	AA117471; C630029K18Rik; Suv4-20h1; Suv420h1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001167885
ORF Size:	2649 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216265).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001167885.1</a> , <a href="#">NP_001161357.1</a>
RefSeq Size:	6131 bp
RefSeq ORF:	2652 bp
Locus ID:	225888
UniProt ID:	<a href="#">Q3U8K7</a>
Cytogenetics:	19 A



[View online »](#)

**Gene Summary:**

Histone methyltransferase that specifically trimethylates 'Lys-20' of histone H4. H4 'Lys-20' trimethylation represents a specific tag for epigenetic transcriptional repression. Mainly functions in pericentric heterochromatin regions, thereby playing a central role in the establishment of constitutive heterochromatin in these regions. KMT5B is targeted to histone H3 via its interaction with RB1 family proteins (RB1, RBL1 and RBL2). Plays a role in myogenesis by regulating the expression of target genes, such as EID3.[UniProtKB/Swiss-Prot Function]