

Product datasheet for **MR216271**

Rnf20 (NM_182999) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rnf20 (NM_182999) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Rnf20
Synonyms:	4833430L21Rik; AW540162; C79397; mKIAA4116
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR216271 representing NM_182999
 Red=Cloning site Blue=ORF Green=Tags(s)


TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTCAGGAATTGAAATAAAAGAGCAGCTGGAGAGCCGGGCACATCTATGCCTCCAGAGAAGAAGACAG
 CTGTTGAAGATTCAGGGACCACCGTGGAAACAATTAAGTTAGGAGGTGTCTTCAACGGAGGAATTAGA
 CATTGCAACTGCAATCCAAGAATCGTAACTAGCAGAAATGCTAGATCAGCGGACGGCCATTGAAGAT
 GAGCTTCGGGAACACATTGAAAAGCTGGAGCGACGCCAGGCCACTGATGATGCCTCACTGTTGATTGTGA
 ACCGGTACTGGAGCCAGTTTGTGAAAACATCCGCATCATCCTTAAACGTTATGATCTGGACCAGGGTTT
 GGGGGACCTACTCACAGAAAGGAAAGCCTTGGTTGTGCCAGAACCAGAACCTGACTCTGATAGCAATCAG
 GAACGCAAAGATACCGTGGAGAGGGGATGGCAAGAGCCAGCTTCTCTTCTTCTGCTACTTTGGCTA
 GCAGTTCCAGTGAAGAGATGGAATCTCAGCTTCAGGAGCGTGTAGAGTCTCCCGCGTGTGTGTCCCA
 GATTGTGACTGTCTATGATAAATGCAAGAAAAGTGGATCTTATCCCGGAAACTGAACAGTGGAGAT
 AACCTGATAGTGGAGGAAGCAGTGCAGGAGTTGAACTCCTTCTTGCACAAGAGAATGTGAGGCTTCAGG
 AGCTGACAGACCTCCTTCAGGAGAAGCATCACACCATGTCTCAGGAGTTCTGTAAGTTGCAGGGTAAAGT
 AGAGACAGCTGAGTACGAGTGTCCGTCTTGGAGTCCATGATTGATGACCTGCAGTGGGATATTGACAAA
 ATTCGAAAAGAGAACAGCGTCTCAACCGGCATTTGGCTGAAGTCTAGAGCGGGTGAATTCAAAAGGCT
 ATAAAGTATATGGAGCGGGGAGCAGTCTCTATGGAGGCACAATCACCATCAATGCCCGGAAGTTTGAAGG
 AATGAATGCAGAGCTTGGAGAGAACAAGAGTTGGCTCAGAACCAGCATTGTGAGCTGGAGAAGCTTCGT
 CAAGACTTTGAAGAGTCACTACACAGAATGAAAAGTTGAAGTGAATTCGGGAGTGCAGTGGAGGAAG
 TGGTTAAGGAAACCCAGAGTATCGCTGCATGCAGTACAGTTCTGTCTGTATAACGAGAGCTTGCA
 GTTGAAGCCCACTTGGATGAAGCTCGACCTGCTTTCATGGCACCCGGGTACTACCAGCGCCAGGTT
 GAACTCATTGAGCGAGATGAGGTTAGTCTGCATAGAAGCTGAGAAGTGAAGTCAATTCAGCTGGAGGATA
 CCTTGGCCAGGTTGCAAGGATGAGATGTTGAGGATAGAGTTTGGAGCAGCCCTGCTGCCAATGA
 ACAAGCAGGCCCAATAAACCGGAGATGCGGCATCTCATTAGCAGCCTCCAGAATCATAACCACCAGCTG
 AAGGTTGAAGTCTAAGTATAAACGGAAGCTGAGAGAGGCGCAGTCTGACCTGAACAAGACAGCTCTGC
 GCAGTGGCAGTCCCTCCTGCAGTCTCAGTCTAGTACTGAGGACCCCAAGGACGAGCCACAGAGCTGAA
 GCAGGATCTGAAGACCTAGCCACTCCTCCTCAGCGCTGAAGGCATCTCAGGAGGATGAAGTTAAGTCC
 AAACGGGATGAAGAAGAGAGAGAACGAGAAAGGAGGGAGAAGAAGAGAGCGAGAGCGAGAACGAGAAA
 AGGAAAAGGAGAGAGAACGAGAGAAACAGAAAAGTGAAGAGTCAAGAAAGGAGAGAGATTCTGTTAAGGA
 TAAAGAGAAAGGAAACATGATGATGGAAGGAAGAAGGAAGCAGAAATTATCAAGCAGCTGAAGATTGAA
 CTTAAGAAAGCACAAGAGAGCCAGAAGGAGATGAAACTACTGCTGGATATGTACCCTCTGCCCAAGG
 AGCAGAGAGACAAAGTGCAGCTAATGGCAGCTGAGAAGAAGTCAAAGGCAGAGTTGGAGGATCTCAGACA
 AAGACTCAAGGATCTGGAAGATAAAGGAGAAGAAAGAGAACAAGAAAATGGCTGATGAGGATGCCTTGAGG
 AAGATCCGGGCAGTGGAGGAGCAGATTGAATACCTGCAGAAGAAGCTGGCTATGGCCAAGCAGGAGGAAG
 AGGCTCTCCTGTCTGAGATGGATGTACAGGACAGGCCCTTGGAGACATGCAGGAGCAGAATCCGCTC
 GATGCAGCAGTTGCGGGAGAAGGACGATGCAAACTTCAAGCTCATGTGAGAGCGAATCAAGTCTAATCAG
 ATCCATAAGTTGCTTAAAGAGGAAAAGGAAGTTAGCTGATCAGGTTTTGACACTGAAGACTCAGGTAG
 ATGCCAGTTACAGGTAGTAAGAAAAGTGAAGAGAAGGAGCACCTACTTCAAGTAAACATTGGCACAGG
 GGAGAAGGAGCTGGGCTTAGGACTCAAGCCTTGGAGATGAATAAACGCAAGGCAATGGAGGCGGCCAG
 CTTGCAGATGACCTCAAAGCACAGCTAGAAGTGGCTCAGAAGAAGCTCCATGATTTCCAGGATGAAATCG
 TGGAGAACAGTGTACCAAAGAAAAGGACTTGTCAATTTCAAACGAGCACAGGAAGACATCTTAGGCT
 TCGAAGAAAGCTAGAGACCACCAAGAAGCCAGACAATGTGCCAAGTGTGATGAAATCTGATGGAGGAG
 ATAAAGGACTACAAGGCACGGCTGACCTGCCCTTGTGTAAACATGCGTAAAAGGATGCAGTACTTACCA
 AGTGTTCACGTTTTCTGCTTTGAGTGTGTAAGACGCGCTATGACACTCGCCAGCGCAAGTGTCCCAA
 GTGTAATGCTGCTTTTGGTGCCAATGATTTCCATCGCATCTACATTGGC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_182999.2 , NP_892044.1
RefSeq Size:	4159 bp
RefSeq ORF:	2922 bp
Locus ID:	109331
UniProt ID:	Q5DTM8
Cytogenetics:	4 B1
MW:	113.5 kDa
Gene Summary:	Component of the RNF20/40 E3 ubiquitin-protein ligase complex that mediates monoubiquitination of 'Lys-120' of histone H2B (H2BK120ub1). H2BK120ub1 gives a specific tag for epigenetic transcriptional activation and is also prerequisite for histone H3 'Lys-4' and 'Lys-79' methylation (H3K4me and H3K79me, respectively). It thereby plays a central role in histone code and gene regulation. The RNF20/40 complex forms a H2B ubiquitin ligase complex in cooperation with the E2 enzyme UBE2A or UBE2B; reports about the cooperation with UBE2E1/UBCH are contradictory. Required for transcriptional activation of Hox genes. Recruited to the MDM2 promoter, probably by being recruited by p53/TP53, and thereby acts as a transcriptional coactivator. Mediates the polyubiquitination of PA2G4 leading to its proteasome-mediated degradation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR216271