

Product datasheet for **MR216303L3V**

Cuedc2 (NM_001164294) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cuedc2 (NM_001164294) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cuedc2
Synonyms:	3010002G01Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001164294
ORF Size:	798 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216303).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001164294.1 , NP_001157766.1
RefSeq Size:	1779 bp
RefSeq ORF:	801 bp
Locus ID:	67116
UniProt ID:	Q9CXX9
Cytogenetics:	19 C3
Gene Summary:	Controls PGR and ESR1 protein levels through their targeting for ubiquitination and subsequent proteasomal degradation.[UniProtKB/Swiss-Prot Function]



[View online »](#)