

## Product datasheet for **MR216313LV**

### **Mpzl1 (NM\_001083897) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Mpzl1 (NM_001083897) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mpzl1
Synonyms:	1110007A10Rik; PZR
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001083897
ORF Size:	810 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216313).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001083897.1</a> , <a href="#">NP_001077366.1</a>
RefSeq Size:	2356 bp
RefSeq ORF:	813 bp
Locus ID:	68481
UniProt ID:	<a href="#">Q3TEW6</a>
Cytogenetics:	1 H2.3



[View online »](#)

**Gene Summary:**

Cell surface receptor, which is involved in signal transduction processes. Recruits PTPN11/SHP-2 to the cell membrane and is a putative substrate of PTPN11/SHP-2. Is a major receptor for concanavalin-A (ConA) and is involved in cellular signaling induced by ConA, which probably includes Src family tyrosine-protein kinases. Isoform 2 seems to have a dominant negative role; it blocks tyrosine phosphorylation of MPZL1 induced by ConA. Isoform 1, but not isoform 2, may be involved in regulation of integrin-mediated cell motility (By similarity).[UniProtKB/Swiss-Prot Function]