

Product datasheet for **MR216327L3V**

Dach1 (NM_007826) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dach1 (NM_007826) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dach1
Synonyms:	Dac; Dach
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_007826
ORF Size:	2253 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216327).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007826.2 , NP_031852.1
RefSeq Size:	5225 bp
RefSeq ORF:	2256 bp
Locus ID:	13134
UniProt ID:	Q9QYB2
Cytogenetics:	14 48.25 cM



[View online »](#)

Gene Summary:

Transcription factor that is involved in regulation of organogenesis. Seems to be a regulator of SIX1, SIX6 and probably SIX5. Corepression of precursor cell proliferation in myoblasts by SIX1 is switched to coactivation through recruitment of EYA3 to the SIX1-DACH1 complex. Transcriptional activation seems also to involve association of CREBBP. Seems to act as a corepressor of SIX6 in regulating proliferation by directly repressing cyclin-dependent kinase inhibitors, including the p27Kip1 promoter. Inhibits TGF-beta signaling through interaction with SMAD4 and NCOR1 (By similarity). Binds to chromatin DNA via its DACHbox-N domain. [UniProtKB/Swiss-Prot Function]