

## Product datasheet for **MR216349L4V**

### **Cd226 (NM\_001039149) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Cd226 (NM_001039149) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cd226
Synonyms:	BC051526; DNAM-1; DNAM1; Pta1; TLISA1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001039149
ORF Size:	660 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216349).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001039149.1</a> , <a href="#">NP_001034238.1</a>
RefSeq Size:	2148 bp
RefSeq ORF:	663 bp
Locus ID:	225825
Cytogenetics:	18 E4



[View online »](#)

**Gene Summary:**

Involved in intercellular adhesion, lymphocyte signaling, cytotoxicity and lymphokine secretion mediated by cytotoxic T-lymphocyte (CTL) and NK cell. Cell surface receptor for NECTIN2. Upon ligand binding, stimulates T-cell proliferation and cytokine production, including that of IL2, IL5, IL10, IL13, and IFNG. Competes with PVRIG for NECTIN2-binding. [UniProtKB/Swiss-Prot Function]