

Product datasheet for MR216352

Ifi203 (NM_008328) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ifi203 (NM_008328) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ifi203
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR216352 representing NM_008328 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGATGGCTGAATACAAGAATATTGTTCTTCTAAAAGGATTAGAGAACATGGAGGACTATCAGTTTAGGA
CAGTCAAGTCCTTACTAAGAAAAGAACTAAAACCTACTAAAAAATGCAAGAAGATTATGATAGGATTCA
GCTGGCGGACTGGATTGAGGACAAGTCCCAAGATGCTGGACTAGACAACTGATAAAAGTGTGTGAG
CACATAAAAGATCTTAAAGACCTTGCTAAAAAAGTAAAAACAGAGAAGGCAAAAGTTCAAAAGAAAAAGG
AAGGAAAATGCAAAACGGCTGGGAAAAAAGGGCAAGATGAACTCAGTTCTTCTGAATCTCTTTTCAA
CAATAAAGAATCAAACAAGAGTGTACCGTCTTCAAAGAAAAAGAGAAAACAGATCACAAAACTGAAGGT
GGAAAGAAAAAGAGCTACCCAGGAGCAGGCTCAGCTCCAGAACTTCAGGAACCAACATAAAAAAAG
AGGAAGATTGCCTCCAGAATCCTCATAAGTCTCCACCAACACCATCCAGCAGTTCCTCAAACAAGGCACC
AAGGAGAGGAACTGTACCAAGGAGCCTTCTAGAGAAGAAGGTCACCATCAAGGTCCCAAACAAGTGATG
GTGTTGAAAGTAACAGAACCATTACATATGACTTTGAAGAAACGAAAAGGATGTTTCATGCTACAGTGG
CTACTGAACTGAGTTCCTCAGAGTGAAGGTTTTGACACAGCTCTAATGAGCAAGTTCATCCAGGAAA
GATCATTGCCATACACATTATTTGGTGCAATGGGTTTTCTGGAGATACAGAGCTTCCTGTGTGTC
GATGTGAACATTAACCCAACAATGATTATCTCAAATACACTGAGTGAAAGTGTATTGCAACTCCAAAA
TTTCTTATCTTCTCACAAAGCAAGGGGACATTTGTAATGGAGAATTTGTAGTATTTAAGAAAAGTGA
GAGGCATGAGTGCATTTGCTATGGAATTGGAGATGATACAGGAAAAATGGCAGTGGTGGTTTATGGACGA
CTCACCAATGTCAGGTGTGAACAGGCAGTAACTTAGACTTGTCTGCTTTGAATTGACTTCCACTAAAG
ATGTGTGCCTCCTGAGGTCTGTAAGGCACAGTTACATGCAGGTCAATGAAGGAAAGCCACTCAACCC
AGACTCAGTCAGGAGAACTCTCTGGAACCATACTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >MR216352 representing NM_008328
Red=Cloning site Green=Tags(s)

MMAEYKNIIVLLKGLNEMEDYQFRTVKSLLRKELKLTKKMQEDYDRIQLADWIEDKFPK DAGLDKLIKVCE
 HIKDLKDLAKKLTKEKAVQKKKEGKCKTAGKKKGQDELSSSESLFNNKESNKSVPSSKKKRKQITKTEG
 GKKKKL TQEQAQLPETS GNTNIKKEEDCLQNP HKSPPTPSSSSSNKAPRRGTVPKEPSREEGHHQGPQVM
 VLKVTPEPTYDFEETKRMFHATVATETEF FRVKVFD TALMSKFIPGKIIAISHYIGCNGFLEIYRASCVS
 DVNINPTMIIISNTLSESAIATPKISYLLSQAKGTFVNGEFVVFKKSERHECICYGIGDDTGKMAVVVYGR
 LTNVRCEPGSKLRLVCFELTSTKDVCLLRSVRHSYMQVINNEGKPLNPDSVRRNSLEPYF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_008328

ORF Size: 1227 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_008328.2](#), [NP_032354.2](#)

RefSeq Size: 3619 bp


RefSeq ORF: 1230 bp

Locus ID: 15950

Cytogenetics: 1 80.76 cM

MW: 46.8 kDa

Product images:



Circular map for MR216352