

Product datasheet for **MR216393L3V**

St6galnac6 (NM_001025311) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	St6galnac6 (NM_001025311) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	St6galnac6
Synonyms:	Siat7f; ST6GalNAcVI
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001025311
ORF Size:	939 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216393).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001025311.2 , NP_001020482.1
RefSeq Size:	2288 bp
RefSeq ORF:	942 bp
Locus ID:	50935
UniProt ID:	Q9JIM95
Cytogenetics:	2 B



[View online »](#)

Gene Summary:

Alpha-2,6-sialyltransferase involved in the synthesis of alpha-series gangliosides. Has activity toward GD1a, GT1b and GM1b. Has no activity toward glycoproteins. Responsible for the biosynthesis of DSGG (disialylgalactosylgloboside) from MSGG (monosialylgalactosylgloboside) in kidney.[UniProtKB/Swiss-Prot Function]