

Product datasheet for **MR216394L4V**

PPP4R4 (NM_028980) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ppp4r4 (NM_028980) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ppp4r4
Synonyms:	8430415E04Rik; mKIAA1622
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_028980
ORF Size:	2625 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216394).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_028980.3 , NP_083256.2
RefSeq Size:	3683 bp
RefSeq ORF:	2628 bp
Locus ID:	74521
UniProt ID:	Q8C0Y0
Cytogenetics:	12 E
Gene Summary:	Putative regulatory subunit of serine/threonine-protein phosphatase 4.[UniProtKB/Swiss-Prot Function]



[View online »](#)