

## Product datasheet for **MR216420**

### **Rab11fip4 (NM\_175543) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Rab11fip4 (NM_175543) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Rab11fip4
Synonyms:	A730072L08Rik; mKIAA1821; RAB11-FIP4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>MR216420 representing NM\_175543  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGCGGGTGGCGGGCTGGCGGGCGCCCCGCGGCTCTGCTGCGCTCCGTGCGCCCTGCGCGAAG  
TGTTTCGAGGTGTGCGGCCGCGACCCCGACGGCTTTCTACGCGTGGAGCGGTGCGCGCGCTCGGCCCTGCG  
CTTCGGGCAAGGAGAGGAGGTGGAAAACTTGTGAAGTGTGGATCCCAATGACCTGGGAGAATCAAC  
TTCAAGGATTTTTGCCGCGGAGTATTCGCCATGAAAGGGTGTGAGGAGCTGCTGAAGGACGTGCTGTCGG  
TGGAGAGCGCAGGGACCTGCCATGTTCTCCAGACATCCAGACTGTGTGGAGCAGGGCAGCGACTTCTC  
TGGCTCTACTGATGGCGAGCAACTCCCCAGGAAACCCGACTTCTCCAGGAAGACGAAGAGGAAGCGATG  
ACATTGGCCCTACCTGAGGGACCCAGGAGTTAGACATGGACAGTCCCATGGAGAGCAGCCAAGGCCCGG  
AGGGCTCTGTCAAGGGTGTGGCGAGGAGAAGGAGCCTGAGCTGGGGGGCCTGTTCTGCCAGAGGACAA  
GTCCCTCGTCTCACTCCATCTGTGACCACATCAGACCTTCCACTCACTCCACCGCCTCTTTATCAGC  
AACGAGGAGCAATTTGAAGACTATGGCGAGGGGATGACGTGGACTGTGCTCCCAGCAGCCCATGCCCCG  
ACGATGAGACCAGAACCAATGTCTACTCAGATCTGGGATCTTCAGTGTCTCCAGTGCGGACAGACCCC  
AAGGAAGATGCGACATGCCTACAACAGCGAGCTGCTGGACGTGTACTGCTCTCAGTGTGCAAGAAGATC  
AACCTGCTGAATGACCTGGAGGCCCGGCTGAAGAACCTGAAGGCCAACAGCCCCAACCGGAAGATCTCAA  
GCACAGCCTTTGGACGGCAGCTCATGCACAACAGTAATTCAGCAGCAGTAATGGCAGCACCGAGGACCT  
GTTCCGTGACAGCATCGACTCCTGTGACAATGACATCACAGAGAAGGTGAGCTTCTGGAAAAAAGGTA  
ACAGAACTGGAGAACGACAGCCTGACGAGCGGGGGCCTGAAGAGCAAGCTGAAGCAAGAGAACATGCAGC  
TGGTGCACAGGGTGCACGAGCTGGAGGAGATGGTGAAGGATCAGGAGACCACAGCGGAGCAGGCGCTGGA  
GGAGGAGGCCCGCGACACCGTGAGGTCTACTGCAAGCTGGAGCGGAAAAGAGCACAGAGCTGGAGCTG  
CTCAACACCAGGGTGCAGCAGTTGGAGGAAGAGAACACAGATCTGAGAACGACAGTGGCTCGACTCAAGT  
CACAGACCAGAAAGCTGGACGAGGAGCGGCAGCGCATGTCTGACCGGCTGGAGGACACCAGCCTGCGACT  
GAAGGACGAGATGGACCTTTACAAGCGCATGATGGACAAGCTGCGGCAGAACCGCCTGGAGTTCAGAAAG  
GAGCGGGAGGCGACGCAGGAGCTCATCGAGGACTTGGGAGAGAGCTGGAGCACCTGCAGATGTACAAGC  
TGGACTGCGAGCGGCCAGGTAGGGGCCGAGCTCCTCTGGCCTTGGCGAGTTCAACGCCAGAGCCCCGGA  
GGTGAACCTGAACATGAGGTCAAGCGTCTGAAGCAGGAGAATCACAAGCTGCGGGATCAGAATGATGAT  
TTAAATGGGCAGATCCTGAGCCTCAGCCTCTACGAGGCGAAGAACCTCTTTGCAACCCAGACCAAAGCCC  
AGTCTCTGGCTGCAGAAATCGACACAGCATCTCGAGATGAGCTAATGGAGGCCCTGAAGGAGCAGGAGGA  
GATCAACTTCCGACTGAGACAGTACATGGACAAGATTATCCTGGCTATACTGGACCAACACCCCTCCATC  
TTGGAGATTAACAC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
TGGATTACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR216420 representing NM\_175543  
Red=Cloning site Green=Tags(s)


MAGGAGWAGAPAAALLRSVRRLREVFEVCGRDPDGLRVERVAALGLRFGQGEEVEKLVKCLDPNDLGRIN  
 FKDFCRGVFAMKGCEELLKDVL SVESAGTLPSPDIPDCVEQGSDFSGSTDGEQLPREPDFQEDEEAM  
 TLALPEGPQELDMSPMESSQGPESVKGCGEKEPELGGFLPEDKCLVLTSPVTTSDLSTHSTASLIS  
 NEEQFEDYGEADDVDCAPSSPCPDDETRNTVYSDLGSSVSSSAGQTPRKMRYHAYNSELLDVYCSQCCKKI  
 NLLNDLEARLKNLKANSPNRKISSAFGRQLMHNSNFSSNGSTEDLFRDSIDSCNDITEKVSFLEKKV  
 TELENDSLTSGGLKSKLKQENMQLVHRVHELEEMVKDQETTAEQALEEEARRHREYVYCKLEREKSTEL  
 LNTRVQQLLEEENTDLRTTVARLKSQTEKLEERQMSDRLEDTSRLKDEMDLYKRMMDKLRQNRLEFQK  
 ERATQELIEDLRRELEHLQMYKLDCEPGRGRSSSGLGEFNARAREVELEHEVKRLKQENHKL RDQND  
 LNGQILSLSLYEAKNLFATQTKAQSAAEIDTASRDELMEALKEQEEINFRLRQYMDKIIILAILDHNPSI  
 LEIKH

SGPTRTRPLEQKLI SEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9092\\_h07.zip](https://cdn.origene.com/chromatograms/mm9092_h07.zip)

**Restriction Sites:** SgfI-RsrII

**Cloning Scheme:**



**ACCN:** NM\_175543

**ORF Size:** 1905 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_175543.3](#), [NP\\_780752.1](#)

**RefSeq Size:** 3145 bp

**RefSeq ORF:** 1908 bp

**Locus ID:** 268451


**UniProt ID:** [Q8BQP8](#)

**Cytogenetics:** 11 46.88 cM

**MW:** 71.9 kDa

**Gene Summary:** Acts as a regulator of endocytic traffic by participating in membrane delivery. Required for the abscission step in cytokinesis, possibly by acting as an 'address tag' delivering recycling endosome membranes to the cleavage furrow during late cytokinesis (By similarity). May play a role in differentiation during retinal development, in a Rab11-independent manner. [UniProtKB/Swiss-Prot Function]

### Product images:



Circular map for MR216420