

## Product datasheet for MR216422

### Fbxo44 (NM\_001161852) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Fbxo44 (NM\_001161852) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Fbxo44  
**Synonyms:** 5730411K09; AV001623; FBG3; Fbx6a; FBX30; Fbxo6a  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >MR216422 representing NM\_001161852  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCAGTAGGCAACATCAACGAGCTGCCAGAGAACATCCTGCTGGAAGTGTTCATCCACATCCCAGCTC  
GCCAACTGCTGCTGCGCTGCCGACCGGTCTGCAGCCTCTGGCGAGACCTCATTGACCTGGTTACACTCTG  
GAAGCGCAAGTGCCTTCAAGAGGGCTTCATCACCGAGGACTGGGACCAGCCTGTGGCTGACTGGAAGATC  
TTCTACTTCTGCGCAGCCTCCAGAGGAACCTCCTTACAACCCTTGTGCTGAAGAGGGATTTGAGTTCT  
GGAGCCTGGACGTGAACGGAGGCGATGAATGGAAGGTAGAGGATCTCTCTAAAGACCAGCGAAAGGAATT  
CCCCAATGACCAGGTTTGCAGCCAGGCCAGACTGCGGGTCCAAGTACCAGCTGTGCGTCCAACTCCTGTC  
GTCAGCCCACGCACCACTGGGAACCTTCCAGCCGACCCGGTGATGATTCAGCAGAAAAGCGATGCCAAG  
TGGAGGGAGGTTTCACACACATTCTCCAATATCCTCCCGCGTCCGCTACATCTGGTTTCAGCACGGAG  
GCGTGGACACTCACTACTGGGCCGGTGGTACGGCCCCGAGAGTCACAAACAGCAGCATCACCATCGGGCC  
CCCGCTGCCCTGATGCCTCCTGAGCCCCAGTCCCTCAGAGCCC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


**Protein Sequence:** >MR216422 representing NM\_001161852  
 Red=Cloning site Green=Tags(s)

MAVGNINELPENILLELFIHIPARQLLLRCRPVCSLWRDLIDLVTLWKRKCLQEGFITEDWDQPVADWKI  
 FYFLRSLQRNLLHNPCAEEGFVWFLDVNGDEWKVEDLSKDQRKEFPNDQVCSQARLRVQVPAVRPTPV  
 VSPRTTGNLPAGPGDDSAEKRCQVEGGFTHILQLSSRRPLHLVSARRRGHSLGLRLVRPESHKQQHHHRA  
 PAALMPPEPQSSEP


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001161852

**ORF Size:** 672 bp

|                               |   |
|-------------------------------|---|
| <b>OTI Disclaimer:</b>        | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>  |
| <b>OTI Annotation:</b>        | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.  |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).  |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| <b>RefSeq:</b>                | <a href="#">NM_001161852.2</a> , <a href="#">NP_001155324.1</a>   |
| <b>RefSeq Size:</b>           | 1688 bp   |
| <b>RefSeq ORF:</b>            | 675 bp  |
| <b>Locus ID:</b>              | 230903  |
| <b>UniProt ID:</b>            | <a href="#">Q8BK26</a>  |
| <b>Cytogenetics:</b>          | 4 E2  |
| <b>MW:</b>                    | 26.1 kDa  |
| <b>Gene Summary:</b>          | Substrate-recognition component of the SCF (SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex.[UniProtKB/Swiss-Prot Function]   |