

Product datasheet for **MR216423**

Chpf (NM_001001566) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Chpf (NM_001001566) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Chpf
Synonyms:	1700028N03Rik; AI414328; D1Bwg1363e
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR216423 representing NM_001001566
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGCGGGCGTCGCTGCTGTCCGTGCTGCGGCCGCGGGGCCGTGGCCGTGGGCATCTCTCTGGGCT
 TCACCCTGAGCCTGCTCAGCGTCACCTGGGTGGAGGAGCCTTGCGGACCCGGGCCGCCCAACCCGAGAG
 CTCTGAGCTGCCCGCGCGCAACACCAACGCGCGCGCCGCCCAACTCGGTGCAGCCGGATCCGAG
 CGCGAGAGGCCCGGGGCCGCGCAGGCACCGGTGAGAGCTGGGAGCCTCGTGTCTTGCCCTACCATCCCG
 CGCAGCCAGGCCAGGCCACCAAGAAGGCCGTGAGAACTCGGTATATCAGCACGGAGCTGGGCATCAGGCA
 GAAGCTTCTGGTGGCAGTGTGACCTCACAAAGCCACGTTGCCACACTGGGTGTGGCTGTAAACCGAACT
 CTGGGACACCGACTGGAGCATGTAGTGTCTGACCGGTGCGAGGGGCCCGGACACCTTCAGGCATGG
 CGGTGGTGGCACTGGGCGAAGAGAGGCCATCGGACACCTGCACCTGGCGCTGCGCCACCTGCTGGAGCA
 ACACGGCGATGACTTTGACTGGTTTTCTAGTGCCTGATGCCACCTATACTGAAGCGCATGGACTGGAC
 CGCTAGCTGGCCACCTCAGCCTTGCTTACGAAACCCATCTCTATCTTGGCCGGCCGAGGACTTCATCG
 GTGGAGATACTACCCAGGCCGCTACTGCCACGGGGCTTTGGAGTCTTGCTCTCTCGCACACTGTACA
 GCAACTGCGCCCCACCTGGAAAGCTGCCGCAACGACATCGTCAGTGCTCGCCGGATGAGTGGTTGGGC
 CGCTGCATCCTTGATGCCACAGGCGTGGGCTGTACTGGTGACCACGAGGGAATGCACTACAACTACCTGG
 AACTGAGCCCCGGGGAGCCTGTGCAGGAGGGGGACCCCTCGTTCCGCGAGCGCCTTGACAGCCATCCCGT
 GCGTGACCCTGTGCATGTACAGCTGCACAAAGCTTTTGGCCGCGTGAGCTGGACCGCACGTACCAG
 GAGATTCAAGAATTGCAGTGGGAGATCCAGAATACCAGCCGACTGGCTGCTGATGGGAGAGAGCCTCTG
 CCTGGCCAGTGGCATCCAGCACCGTCTCGCCCTGCCTCAGCCTTTGAGGTTCTGCGCTGGGACTACT
 CACAGAACAATACGCGTTCTCCTGCGCCGATGGCTCTCCCGCTGCCCGTTGCGTGGGGCCGACAGGCT
 GATGTGGCTGACGTCTGGGACAGCCTTAGAGGAGCTCAACCGCGTTACCAGCCAGCGCTGCAGCTCC
 AGAAGCAACAGCTGGTGAACGGTACCAGGCTTTTATCCAGCCGAGGCATGGAGTACACACTAGACCT
 GCAGCTGGAAGCGCTGACACCCAGGTTGGCCGCTGGCCCTCACCCGAGGTTGCAGCTCCTTCGGCCC
 TTGAGCCGAGTGGAGATCTTGCTGTACCCTATGTCACCGAGGCTTCTCGGCTCACTGTGCTACTGCCGC
 TGGCTGCAGCGGAACGAGACCTGGCTTCTGGCTTCTTAGAAGCCTTTGCCACTGCAGCCCTGGAACCTGG
 TGATGCAGCAGCCTTGACCCTGCTGCTGTATGAGCCACGCCAGGCCAGCGGGCAGCCCACTCAGAC
 GTCTTCGCACCTGTCAAGGCCACGTGGCAGAGCTAGAGCGGCGTTCCCTGGTGCCCGGTGCCCTGGC
 TCAGTGTGCAGACAGCAGCGCCCTCCTCACTGCGTCTCATGGATCTGCTGTCCAAGAAGCACCCACTAGA
 CACTCTGTTCTGCTGGCCGGCCAGACACGGTACTCACACCTGATTTCCCTGAACCGCTGCCGCATGCAT
 GCCATCTCTGGCTGGCAGGCTTCTTCCCATGCACCTCCAGGCTTCCACCCTGCTGTGGCTCCACCTC
 AGGGCCCTGGGCCACCAGAGCTGGGCCGTGACACCGGTCACTTTGATCGCCAGGCTGCCAGTGAGGCATG
 CTTCTACAACCTCCGACTATGTGGCGGCCGTGGCCGGCTGGTGGCGGCCCTCGGAGCAGGAGGAGGAGCTG
 CTGGAGAGCCTGGATGTGTACGAGTTGTTCTGCGCTTCTCAACTTGACAGTGTGAGAGCAGTAGAGC
 CAGCCTTGCTGCAGCGCTACCGGGCCAGCCGTGACGTGCACGGCTCAGTGAAGACCTTTACCACCGCTG
 TCGCCAGAGCGTACTTGAGGGCTTGGTCCCGCACCCAGCTTGCCATGCTGCTCTTTGAGCAGGAACAG
 GGAACAGCACC

AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR216423 representing NM_001001566
Red=Cloning site Green=Tags(s)


MRASLLLSVLRPAGPVAVGISLGFTLSLLSVTWVEEPCGPGPPQPGDSELPARGNTNAARRPNSVQPGSE
RERPGAGAGTGESWEPRVLPYHPAQPGQATKKAVRTRYISTELGIRKLLVAVLTSQATLPTLGVAVNRT
LGHRLEHVVFLTGARGRTPSGMAVVALGEERPIGHHLALRHLLLEQHGDFFDWFFLVPDATYTEAHGLD
RLAGHLASATHLYLGRPQDFIGGDTTPGRYCHGGFGVLLSRTLLQQLRPHLESCRNDIVSARPDEWLG
RCILDATGVGCTGDHEGMHYNLELSPGEPVQEGDPRFRSALTAHPVRDPVHMYQLHKAFARAELDRYQ
EIQELQWEIQNTSRLAADGERASAWPVGIPAPSRPASRFEVLRWDYFTEQYAFSCADGSPRCPLRGADQA
DVADVLGTALEELNRRYQPALQLKQQLVNGYRRFDPARGMEYTLDLQLEALTPQGGRWPLTRRVQLLRP
LSRVEILPVPYVTEASRLTVLLPLAAAERDLASGFLEAFATAALEPGDAAALTTTTLLEPRQAQRAAHS
VFAPVKAHVAELERRFPGARVPWLSVQTAAPSPLRLMDLLSKKHPLDTLFLLAGPDTVLTDFLNRCRMH
AISGWQAFPMHFQAFHPAVAPPQGGPPELGRDTGHFDRQAASEACFYNSDYVAARGRLVAASEQEEEL
LESLDVYELFLRFNLHVLRAVEPALLQRYRAQPC SARLSEDLYHRCRQSVLEGLGSRTQLAMLLFEQEQ
GNST

SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9102_e07.zip

Restriction Sites: Sgfl-RsrII

Cloning Scheme:



ACCN: NM_001001566

ORF Size: 2322 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001001566.3](#), [NP_001001566.1](#)

RefSeq Size: 3107 bp

RefSeq ORF: 2325 bp

Locus ID: 74241


UniProt ID: [Q6IQX7](#)

Cytogenetics: 1 39.14 cM

MW: 85.5 kDa

Gene Summary: Has both beta-1,3-glucuronic acid and beta-1,4-N-acetylgalactosamine transferase activity. Transfers glucuronic acid (GlcUA) from UDP-GlcUA and N-acetylgalactosamine (GalNAc) from UDP-GalNAc to the non-reducing end of the elongating chondroitin polymer (By similarity). Isoform 2 may facilitate PRKN transport into the mitochondria. In collaboration with PRKN, isoform 2 may enhance cell viability and protect cells from oxidative stress (By similarity). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR216423