

Product datasheet for MR216430

Dbnnd1 (NM_028146) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Dbnnd1 (NM_028146) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Dbnnd1
Synonyms: 2810427I04Rik; D8Ertd590e
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR216430 representing NM_028146
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCGCTGTGCCGAAGCCGGACTCAGTAGTCACCAGAGCAGTTGGGGCACCCGGGCTGTCGCCGCCG
TGCCACACCCCGCGGACCCATTGGCCGCATGGAGTCTCCCGAGGGCGTGGTCCGGGAGAGATCGT
TAAGGACGTCAAGGTACCACAGGCTGCCCTGAATGTTTCAGCTCATGAGACAGGGGACATGTGCCGTACA
CCTGTGGCAGAGGAGGAGGAGGTTGGCATCCAATACCAGCACCAGGGTTCCTGCAGGTCACAGAGA
GGAGGCAGCCTCTGAGCAGCGTCTCTCTCTGGAGGTCCACTTCGATCTCTCTGGACCTCACTGAGCTGAC
GGACATGTCTGACCAGGAGTTGGCCGAGGTCTTTGCGGACTCCGATGACGAAAACCTTGCCACTGAGTCA
CCAGCAGGTCTGCACCCACTATCTCGAGCCAGCTGCCTACGTTACCATCCTGGACGAAAACAAGGGCGG
AGCAGAACCGTGAGAAACAGCCCCCAGTGACCCAGAGCGGCAAGGGACAATTGTGGACACATTCCTCAC
TGTGGAGGAGCCAAGGAAGAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR216430 representing NM_028146
 Red=Cloning site Green=Tags(s)

MRCARSPDSVVTRAVGAPGLSPPLTPPADPIGRMESPEGAGPGEIVKDVKVPQAALNNSAHETGDMCRT
PVAEEEEVGIPIAPGFLOVTEREQPLSSVSSLEVHFDLLDL TELTDMSDQELAEVFASSDENLATES
PAGLHPLSRASCLRSPSWTKTRAEQNREKQPPSDPERQGTIVDTFLTVEEPKED

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Chromatograms: https://cdn.origene.com/chromatograms/mm9078_h10.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_028146

ORF Size: 582 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_028146.4](#), [NP_082422.3](#)


RefSeq Size: 1681 bp

RefSeq ORF: 585 bp

Locus ID: 72185

UniProt ID: Q9CZ00
 Cytogenetics: 8 72.1 cM
 MW: 20.9 kDa

Product images:



Circular map for MR216430