

Product datasheet for **MR216437**

Slc27a3 (NM_011988) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Slc27a3 (NM_011988) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Slc27a3
Synonyms:	Acsvl3; FATP-3; FATP3; Vlcs-3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR216437 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGGCCCTCTGCTGCTGCTGCCCTTACTCTTGCTCCTGCCTCTGCTGCTGAAGCTAGACGTCTGGC
CTCAGCTGCGCTGGCTGCCGGCCGACCTGGCCTTCACAGTTCGCGCCCTCCGCTGCAAAAAGGGCCCTTCG
AGCACGCGCGCTGGCCGAGCGGCCGCCGACCCGAAAGCTCTGAGAGCGGGTGCACTCTGGCCTGGCGT
CTCGCTACCTGGCCCGGAGCAGCCGACACACACCTTCTCATCCACGGCGCGCAGCGCTTTAGCTACG
CGGAGGCTGAGCGGAGAGCAACCGGATTGCTCGCGCTTTCTGCGCGCACGGGGCTGGACCGGGGGCCG
CCGAGGCTCGGGCAGGGCAGCACTGAGGAAGGCGCACGCGTGGCGCTCCGGCTGGAGATGCGGCTGCT
GGAGGGACGACCGCCCTCTGGCACCCGGGGCACCGTGGCGCTGCTCTCCAGCGGGCCCGGATT
TCCTTTGGATTTGGTTCCGACTGGCCAAAGCTGGCCTGCGCACGGCTTTGTGCCACCGCTTTACGCCG
AGGACCCCTGCTGCACTGCCTCCGAGCTGCGGTGCGAGTGGCTCGTCTGGCCACAGAGTTCTGGAG
TCCTTGAGCCGACCTGCCGGCTTGAGAGCCATGGGGCTCCACCTATGGGCGACGGGCCCTGAAACTA
ATGTAGCTGGAATCAGCAATTTGCTATCGGAAGCAGCAGACCAAGTGGATGAGCCAGTGCCGGGGTACCT
CTTGCCCCCAGAACAATGGACACCTGCCTGTACATCTTACCTCTGGCACTACTGGCTGCCAAG
GCTGCTCGAATCAGTCATCTGAAGTTCTACAGTGCCAGGGATTCTACCATCTGTGTGGAGTCCACCAGG
AGGACGTGATCTACCTCGCACTCCCACTGTACCACATGTCTGGCTCCCTTCTGGGCATTGTGGGTGCTT
GGGCATTGGGGCCACCGTGGTGTGAAACCAAGTTCTCAGCTAGCCAGTTCTGGGACGATTGCCAGAAA
CACAGGTGACAGTGTCCAGTACATTTGGGAGTTGTGCCGATACCTCGTCAACCAGCCCCGAGCAAGG
CAGAGTGTGACCATAAGGTGCGCTTGGCAGTGGGAGTGGTTGCGCCAGACACCTGGGAGCGTTTCT
GCGGCGATTTGGACCTCTGCAGATACTGGAGACGTATGGCATGACAGAGGGCAACGTAGCTACGTTCAAT
TACACAGGACGGCAGGGTGCAGTGGGGCAGCTTCTGGCTTTACAAGCACATCTTCCCCTTCTCCTTGA
TTCGATACGATGTCATGACAGGGGAGCCTATTCGGAATGCCAGGGGCACTGCATGACCACATCTCCAGG
TGAGCCAGGCCTACTGGTGGCCCCAGTGAGCCAGCAGTCCCCCTTCTGGGCTATGCTGGGGCTCCGGAG
CTGGCCAAGGACAAGCTGCTGAAGGATGTTCTGGTCTGGGACGTTTTCTTCAATACTGGGGACCTCT
TGGTCTGTGATGAGCAAGGCTTCTTCACTTCCACGATCGTACTGGAGACACCTCAGGTGGAAGGGAGA
GAATGTGGCCACAACGAAGTGGCTGAGGTCTTGGAGACCCTGGACTTCTTCAAGGAGTGAACATCTAT
GGAGTACCGGTGCCAGGGCACGAAGGCAGGGCAGGCATGGCGGCTTGGCTCTGGGCCCCCGCAGGCTC
TGAACCTGGTGCAGCTCTACAGCCATGTTTCTGAGAACTTGCCACCGTATGCCCGACCTCGGTTTCTCAG
GCTCCAGGAATCTTTGGCCACTACTGAGACCTTCAAACAGCAGAAGGTTAGGATGGCCAATGAGGGCTTT
GACCCAGTGTACTGTCTGACCCACTCTATGTTCTGGACCAAGATATAGGGGCCTACCTGCCCTCACAC
CTGCCCGTACAGTGCCCTCCTGTCTGGAGACCTTCAATC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR216437 protein sequence
Red=Cloning site Green=Tags(s)


MAALLLLPLLLLLPLLLKLDVWPQLRWLPADLAFTVRALRCKRALRARALAAAAADPESSESGCSLAWR
 LAYLAREQPTHFLIHGAQRFSYAEAERESNRIARAFLRARGWTGRRGSGRGSTEEGARVAPPAGDAAA
 GGTTAPPLAPGATVALLLPAGPDFLWIWFLAKAGLRTAFVPTALRRGPLLHCLRSCGASALVLATEFLE
 SLEPDLPALRAMGLHLWATGPETNVAGISNLLSEAADQVDEPVPGYLSAPQNMIDTCLYIFTSGTTGLPK
 AARISHLKVLQCQGFYHLCGVHQEDVIYLALPLYHMSGSLLGIVGCLGIGATVVLKPKFSASQFWDDCQK
 HRVTVFYIGELCRYLVNQPPSKAECDHKVR LAVGSLRPDPTWERFLRRFGPLQILETYGMTEGNVATFN
 YTGRQAVGRASWLKXIFPFSLIRYDVTGEPIRNAQGHCM TSPGEPGLLVAPVSQQSPFLGYAGAPE
 LAKDKLLKDVFWSGDVFNTGDLLVCDEQGF LHFHDRTGDTFRWKGENVATTEVAEVLETDFLQEVNIY
 GVTVPGHGGRAGMAALALRPPQALNLVQLYSHVSENLPYARPRFLRLQESLATTETFKQKVRMANEGF
 DPSVLS DPLYVLDQDIGAYLPLTPARYSALLSGDLRI

SGPTRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-RsrII

Cloning Scheme:



ACCN: NM_011988

ORF Size: 2004 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_011988.1](#), [NM_011988.2](#), [NM_011988.3](#), [NP_036118.2](#)

RefSeq Size: 2284 bp

RefSeq ORF: 2004 bp

Locus ID: 26568


UniProt ID: [O88561](#)

Cytogenetics: 3 F1

MW: 72.9 kDa

Gene Summary: Has acyl-CoA ligase activity for long-chain and very-long-chain fatty acids. Does not exhibit fatty acid transport activity.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR216437