

Product datasheet for **MR216437L3V**

Slc27a3 (NM_011988) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc27a3 (NM_011988) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slc27a3
Synonyms:	Acsvl3; FATP-3; FATP3; Vlcs-3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_011988
ORF Size:	2004 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216437).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011988.2 , NP_036118.2
RefSeq Size:	2284 bp
RefSeq ORF:	2004 bp
Locus ID:	26568
UniProt ID:	O88561
Cytogenetics:	3 F1
Gene Summary:	Has acyl-CoA ligase activity for long-chain and very-long-chain fatty acids. Does not exhibit fatty acid transport activity.[UniProtKB/Swiss-Prot Function]



[View online »](#)