

Product datasheet for **MR216447**

Kmt5c (NM_001115018) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Kmt5c (NM_001115018) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Kmt5c
Synonyms: BC024816; Suv4-20h2; Suv420h2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR216447 representing NM_001115018
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGACCAATGGGGCCAAGATTGTGTCTACCCGTGCTTGAAGAAGAATGAGAAGCTGGAGCTGCTGG
TGGGCTGCATTGCAGAGCTGCGTGAAGAGGATGAAGACCTGCTTAGGGCCGGTGAAGACGACTTCAGCAT
CATGTACTCCACCGAAAAGGAGTGCCAGCTGTGGCTTGGCCAGCTGCCTTCATCAACCATGACTGC
AAACCCAACTGCAAGTTTGTGCCCTCAGATGGGAATACTGCCTGTGTGAAGGTGCTACGGGACATTGAGC
CAGGGGATGAAGTGACTTGTCTATGGTGAGGGCTTCTCGGGCAGAAGAATGAGCACTGTGAATGCTA
TACCTGTGAGAGGAAAGGCGAGGGAGCTTTCAGACTTCAGCCAGGGAGCCTGAGCTGCGACCAAAGCCC
CTAGACAAGTATGAGCTCCGAGAGACAAAGAGGCGTCTGCAGCAGGGCTTGGTCAGCAGCCAGCAGAGCC
TGATGAGCCGGTGGCCCTGCTCCACCTGTCCCACTGCGCCCTGACCCATTCTGTGCTGCCTGCCAGCC
CTCATGCCTGCTGCCTGCCAGCCCCACATGGATTACTTGCCCTTTGGCTCCAGCGGGCGCCTCAACCC
CAGCCCATAGTACCTCCAGGAAGCGCCACCGCCGCCGCCAGGATCCGCCAGGCTTCATTGCCTC
CTGTCTCCGCACTGCCTGTGTCCCTGCACCGATGGGGAGGCTGTGGTCCCACTGCCAGCTACGTGC
TGAGGCCATGGTGACCCTACACCTGCGCCCCAGACTCGCTGGACCCACAACAGGACTGGTACTGGGCC
CGGCGCTATGGGCTACCTTCTGTGGGGCGTGTGGAACCTACCCGCTAGCCCCAGCCCTGCCAGCTGCTC
CTGCCCTGCTGGGAACCCAGGTCCTGTCCCACTCCTGATTTTCATCCCTAAGCAGGCCCTTGCTTTGC
TCCGTTTGCCACCTAAGCGCCTCCGGCTAGTGGTCAGTCATGGTTCTATTGATCTGGACATCAATAGT
GGTGAGCCA

ACGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online >](#)

Protein Sequence: >MR216447 representing NM_001115018
 Red=Cloning site Green=Tags(s)


METNGAKIVSTRAWKKNEKLELLVGCIAELREEDDLLRAGENDFSIMYSTRKRSACLWLGPAAFINHDC
 KPNCKFVPSDGNACVKVLRDIEPGDEVTCFYGEGFFGEKNEHCCEYTCERKGEAFRLQPREPELRPKP
 LDKYELRETKRRLQQGLVSSQQSLMSRWACSHLSPLRPDPFCAACQPSCLLPASPHMDYLLPLWLQRAPQP
 QPIVPPRKRHRRRRPRIRQASLPPVLRACVPLHRWGGCGPHCQLRAEAMVTLHLRPQTRWTPQQDWYWA
 RRYGLPSVGRVELTRLAPALPAAPAPAGNPGVPVTPDFIPKQALAFAPFCPPKRLRLVSHGSIDLIDINS
 GEP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001115018

ORF Size:	1404 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq Size:	2191 bp
RefSeq ORF:	1407 bp
Locus ID:	232811
UniProt ID:	Q6Q783
Cytogenetics:	7 A1
MW:	53.6 kDa
Gene Summary:	Histone methyltransferase that specifically trimethylates 'Lys-20' of histone H4. H4 'Lys-20' trimethylation represents a specific tag for epigenetic transcriptional repression. Mainly functions in pericentric heterochromatin regions, thereby playing a central role in the establishment of constitutive heterochromatin in these regions. KMT5C is targeted to histone H3 via its interaction with RB1 family proteins (RB1, RBL1 and RBL2).[UniProtKB/Swiss-Prot Function]