

Product datasheet for **MR216454L4V**

Wdr82 (NM_029896) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Wdr82 (NM_029896) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Wdr82
Synonyms:	9430077D24Rik; CDW5/WDR82
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_029896
ORF Size:	942 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216454).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_029896.1 , NP_084172.1
RefSeq Size:	4313 bp
RefSeq ORF:	942 bp
Locus ID:	77305
UniProt ID:	Q8BFQ4
Cytogenetics:	9 F1



[View online »](#)

Gene Summary:

Regulatory component of the SET1 complex implicated in the tethering of this complex to transcriptional start sites of active genes. Facilitates histone H3 'Lys-4' methylation via recruitment of the SETD1A or SETD1B to the 'Ser-5' phosphorylated C-terminal domain (CTD) of RNA polymerase II large subunit (POLR2A). Component of PTW/PP1 phosphatase complex, which plays a role in the control of chromatin structure and cell cycle progression during the transition from mitosis into interphase. Possible role in telomere length maintenance and in mRNA processing (By similarity).[UniProtKB/Swiss-Prot Function]