

## Product datasheet for **MR216458L4V**

### Zrsr1 (NM\_011663) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Zrsr1 (NM_011663) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Zrsr1
Synonyms:	D11Ncvs75; Irlgs2; SP2; U2af1-rs1; U2afbp-rs
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_011663
ORF Size:	1284 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216458).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_011663.3</a> , <a href="#">NP_035793.1</a>
RefSeq Size:	4492 bp
RefSeq ORF:	1287 bp
Locus ID:	22183
UniProt ID:	<a href="#">Q64707</a>
Cytogenetics:	11 14.24 cM
Gene Summary:	Plays a role in splicing of the U12-type introns (PubMed:29617656). Implicated also in removal of U2 introns positioned adjacent to a U12 intron (PubMed:29617656). [UniProtKB/Swiss-Prot Function]



[View online »](#)