

Product datasheet for **MR216477L4V**

KIhI38 (NM_177755) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	KIhI38
Synonyms:	8230402K04Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_177755
ORF Size:	1743 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR216477).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_177755.3 , NP_808423.1
RefSeq Size:	2212 bp
RefSeq ORF:	1746 bp
Locus ID:	268807
UniProt ID:	Q8BSF5
Cytogenetics:	15 D1

