

Product datasheet for **MR216508**

Pgm3 (NM_001163746) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pgm3 (NM_001163746) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Pgm3
Synonyms:	2810473H05Rik; Agm1; BB187688; C77933; PAGM; Pgm-3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR216508 representing NM_001163746
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGATCTGGAGGCTGTGTGTAAGCGCTCGGCACTCCATGCCAAGCCCCAGGGGCTCATTCTTCAGTATG
 GGACTGCTGGATTTCGGACCAATGCCAGCATCTCGATCATATCATGTTTCGAATGGGATTGCTAGCTGT
 CCTGCGGTCAAAGCAAACAAGATCTACCATTGGAGTCATGGTTACGGCCTCACATAACCCCTGAGGAAGAC
 AATGGTGTGAAATTGGTTGATCCTTTGGGGGAAATGTTGGCCCTTCTGGGAGGAGCATGCCACCTGTC
 TGGCCAGTGGGAGGAGCAGGACGTGGGCAAGTCTTGTGCCATCGTGGAGAAAGAAGCTGTGGATCT
 GACACAGACCGCCTTCGTAGTCATTGCCAGAGACACCAGGCCAGCAGTGAGAACTTTCACAATCTGTG
 ATAGATGGCGTGACTGTTTTAGGAGTCAATCCATGATTACGGCCTGTTGACGACACCTCAGTTGCACT
 ACATGGTGTATTGCCGAATCCGGAGGCCGTATGGACAGGCCACCGTAGAAGGCTACTGCCAGAAGCT
 CTCCAAGCTTTCTAGACCTTACGAATCAGGTTTCTGCAGTGGAGATGTAAGAGGTCAGTTAAGGTT
 GACTGTGCCAACGGCATAGGGCCCTTGAAGCTAAGAGAAATGGAACACTACTTCTCCCGGGCCTGTCCG
 TTCTGCTGTTAATGATGGGACCAAGGGCGGCTCAATCACCTGTGCGGTGCTGACTTTGTGAAAAGTCA
 ACAGAAACCCCCACAGGAATTGAAATGAAGTCCGGTGAGAGATGCTGCTCCTTCGATGGGGATGCGGAC
 AGGATCGTGTATTACTACTGTGATGCAGATGGTCACTTTCATCTCATAGATGGAGACAAGATAGCGACGT
 TAATTAGCAGTTTCCTTAAAGAGCTACTGTTGGAGATTGGAGAAAGTGTGAACCTCGGAGTCGTACAGAC
 AGCATATGCAAAATGGAAGTTCACGCGGTACCTTGAAGAAGTTATGAAGGTACCTGTGTATTGCACAAA
 ACTGGTGTAAACATTTGCATCACAAGGCTCAAGAGTTTGACATCGGAGTTTATTTTGAAGCGAACGGGC
 ATGGAACAGCACTGTTTCAGTGAAGCAGTTGAAGTGAAGTAAAAAGACTAGCCCAAGAAGTACAGCATGG
 AAAAGGAAAAGCAGCCAGGACGCTTGCAGCATCATCGACTTATTCAACAGGTCGCGGACAGGAGATT
 ATTAGCACCGATGCTGAGAGACAAGCAGTCACTCCACCAGGACTCCAAGAGGCAATCAATGACCTGG
 TGAAGAAATACACACTTGCCCGAGCTTTGTACGGCCCTCTGGTACTGAAGACATTGTCCGAGTATATGC
 AGAAGCCAACCTACAAGAAAGCGCAGACAGGCTTGCCATGAAGTGAAGTTGCTGGTGTCCAGCTGGCT
 GGAGGCATTGGAGAAAGACCACAACCAACTTTC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR216508 representing NM_001163746
 Red=Cloning site Green=Tags(s)

MDLEAVCKRSALHAKPQGLILQYGTAGFRNTAQHLDHIMFRMGLLAVLRSKQTRSTIGVMVTASHNPEED
 NGVKLVDPLEMLAPSWEEHATCLASAEQDVRQVLAATIVEKEAVDLTQTAFVVIARDTRPSSEKLSQSV
 IDGVTVLGGQFHGYLLTTPQLHYMVYCRNSGGRYGQATVEGYCQKLSKAFVDLTNQVSCSGDVKRSVKV
 DCANGIGALKLREMEHYFSRGLSVLLFNDGTQGRNLHLCGADFVKSQKPPQGGIEMKSGERCCSFDGDAD
 RIVYYYCDADGHFHLIDGDKIATLISSFLKELLEIGESVNLGVVQTAYANGSSTRYLEEVMKVPVYCTK
 TGVKHLHHKAQEFDIGVYFEANGHTALFSEAVEVKIKRLAQELDDGKGAARTLASIIDLFNQVADRRV
 ISTTDAERQAVTTPGLQEAINDLVKKYTLARAFVRPSGTEDIVRVYAEANSQESADRLAYEVSLLVFQLA
 GGIGERPQPTF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:



ACCN: NM_001163746

ORF Size: 1503 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001163746.1](#), [NP_001157218.1](#)

RefSeq Size: 4459 bp

RefSeq ORF: 1506 bp

Locus ID: 109785


UniProt ID: [Q9CYR6](#)

Cytogenetics: 9 46.58 cM

MW: 55.4 kDa

Gene Summary: Catalyzes the conversion of GlcNAc-6-P into GlcNAc-1-P during the synthesis of uridine diphosphate/UDP-GlcNAc, a sugar nucleotide critical to multiple glycosylation pathways including protein N- and O-glycosylation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR216508