

Product datasheet for MR216511L4

Nmur1 (NM_010341) Mouse Tagged Lenti ORF Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | Nmur1 (NM_010341) Mouse Tagged Lenti ORF Clone |
| Tag: | mGFP |
| Symbol: | Nmur1 |
| Synonyms: | FM-3; Gpr66; NmU-R; NMU-R1; NMU1R |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| E. coli Selection: | Chloramphenicol (34 ug/mL) |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR216511). |
| Restriction Sites: | SgfI-MluI |
| Cloning Scheme: | |

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.



[View online »](#)

Locus ID: 14767
UniProt ID: [O55040](#)
Cytogenetics: 1 C5
Gene Summary: Receptor for the neuromedin-U and neuromedin-S neuropeptides.[UniProtKB/Swiss-Prot Function]