

Product datasheet for **MR216520L3V**

Bri3bp (NM_029752) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Bri3bp (NM_029752) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Bri3bp
Synonyms:	2410150118Rik; A1841257; AW742481
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_029752
ORF Size:	759 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216520).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_029752.2 , NP_084028.1
RefSeq Size:	7001 bp
RefSeq ORF:	762 bp
Locus ID:	76809
UniProt ID:	Q8BXV2
Cytogenetics:	5 G1.1
Gene Summary:	Involved in tumorigenesis and may function by stabilizing p53/TP53.[UniProtKB/Swiss-Prot Function]



[View online »](#)