

Product datasheet for MR216533

Dpy30 (NM_001146222) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Dpy30 (NM_001146222) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Dpy30
Synonyms: 2810410M20Rik; C87842
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR216533 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGGATCGCC

ATGGAGTCGGAGCAGATGCTGGAGGGACAGACACAGGTTGCAGAAAACCCTCACTCTGAGTACGGACTCA
 CAGACAGCGTTGAGAGAATAGTCGAAAATGAGAAGATTAATGCAGAAAAGTCATCGAAACAGAAGGTGGA
 TCTACAGTCCTTGCCACCCGTCCTACCTGGACCAGACAGTTGTGCCTATCTTATTACAGGGACTTGCT
 GTGCTTGCCAAGGAAAGACCACCAAATCCCATCGAGTTCCTAGCATCTTATCTTTAAAAACAAGGCGC
 AGTTTGAAGATCGAAAT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR216533 protein sequence
 Red=Cloning site Green=Tags(s)
 MESEQMLEGQTQVAENPHSEYGLTDSVERIVENEKINAEEKSKQKVDLQSLPTRYLDQTVPILLQGLA
 VLAKERPPNPIEFILASYLLKNAQFEDRN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001146222

ORF Size: 300 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001146222.1](#), [NP_001139694.1](#)


RefSeq Size: 790 bp

RefSeq ORF: 300 bp

Locus ID: 66310
UniProt ID: [Q99LTO](#)
Cytogenetics: 17 E2
MW: 11.2 kDa

Gene Summary: As part of the MLL1/MLL complex, involved in the methylation of histone H3 at 'Lys-4', particularly trimethylation. Histone H3 'Lys-4' methylation represents a specific tag for epigenetic transcriptional activation. May play some role in histone H3 acetylation. In embryonic stem (ES) cells, plays a crucial role in the differentiation potential, particularly along the neural lineage, regulating gene induction and histone H3 'Lys-4' methylation at key developmental loci, including that mediated by retinoic acid. Does not affect ES cell self-renewal. May also play an indirect or direct role in endosomal transport.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR216533