

Product datasheet for **MR216558**

Txn14a (NM_001042408) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Txn14a (NM_001042408) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Txn14a
Synonyms:	D18Wsu98e; Dim1; ENSMUSG00000057130; Txn14; U5-15kD; U5-15kDa
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR216558 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s) TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC GCC CGATCGCC ATGTCGTACATGCTTCCGCATCTGCACAATGGCTGGCAGGTAGACCAGGCCATCCTTTCGGAGGAGGACC GCGTGGTCGTCATTCTGTTTCGGACACGACTGGGACCCCACTTGCATGAAGATGGACGAGGTTCTGTACAG CATCGCCGAAAAGAAATGGAAGATAGTGGGAGATCTTCTCATCTTGT ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT ACAAGGATGACGACGATAAGGTTTAA
Protein Sequence:	>MR216558 protein sequence Red=Cloning site Green=Tags(s) MSYMLPHLHNGWQVDQAILSEEDRVVIRFGHDWDPTCMKMDEVLYSIAEKKWKIVGDLPHLV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Restriction Sites:	Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001042408

ORF Size: 192 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001042408.2](#)

RefSeq Size: 1659 bp

RefSeq ORF: 192 bp
 Locus ID: 27366
 Cytogenetics: 18 53.34 cM
 MW: 7.4 kDa
 Gene Summary: Plays role in pre-mRNA splicing as component of the U5 snRNP and U4/U6-U5 tri-snRNP complexes that are involved in spliceosome assembly, and as component of the precatalytic spliceosome (spliceosome B complex).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR216558