

Product datasheet for **MR216610L3V**

F2r12 (NM_010170) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	F2r12 (NM_010170) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	F2r12
Synonyms:	F730031A08; PAR3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_010170
ORF Size:	1107 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216610).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010170.4 , NP_034300.3
RefSeq Size:	2464 bp
RefSeq ORF:	1110 bp
Locus ID:	14064
UniProt ID:	O08675
Cytogenetics:	13 50.38 cM
Gene Summary:	High affinity receptor for activated thrombin coupled to G proteins that stimulate phosphoinositide hydrolysis. May play a role in platelets activation.[UniProtKB/Swiss-Prot Function]



[View online »](#)