

Product datasheet for **MR216611**

Arhgap25 (NM_001037727) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Arhgap25 (NM_001037727) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Arhgap25
Synonyms:	A130039I20Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR216611 representing NM_001037727
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTCTCTCAAGCTGCCAAGGAAGTGGGATTTCAATTTGAAAGCGGAGGCTTCAAAAATAGCTCGGTGCGA
 GGAGTGTGATGACCGGGGAGCAGATGGCGGCCTTTCACCCACCCACCACTCCCAACCCCTGAAAAGGCC
 AATTAAGGTGGGCTGGCTGAAGAAGCAGAGTCCATCGTAAAAACTGGCAGCAGCGTACTTTGTGTTG
 AGAGCTCAGCAACTCTACTACTACAAGGACGAAGAGGACTCAAAGCCACAGGGCTGCATGTATCTGCCAG
 GCAGTACAGTCAAAGAAATTGCCACAAACCCGAAGAAGCTGGGAAATTTGTCTTTGAAGTTATTCCAGC
 CTCAAGTGACCAGAACCATTGGACAAGACTCCTACGTCCTCATGGCCAGCTCCAGGTAGAGATGGAG
 GAGTGGTTAAGTTCCTCAGGAGAGTGTGTCACGCCCTCTGGAGCGGTGTTTGGCCAGCGTCTGGATG
 AGACTGTGGCCTATGAGCAGAAGTTTGGCCCTCACCTGGTGCCCATCTTGGTGGAGAAGTGCCTGAGTT
 CATCCTGGAACATGGTGTGAGTGAAGAGGGCATCTTCCGCTGCCCGGGCAGGACAACCTTGTGAAGCAG
 CTGAGAGATGCTTTCGATGCGGGGAGCGGCCCTCTTTGACAGGGACACAGATGTGCACACGGTGGCAT
 CTCTACTAAAGCTCTACCTCCGAGACTGCCAGAGCCTGTGGTTCCTTGAGTCAAGTATGAAGGGTCTCT
 GCTCTGTGGGAGCTCATGAACGCAGATGAGGCAAAGGCTCAGCAGGAGTTGGTGAAGCAGCTTTCTACC
 CTTCCCGAGACAACATAACCTCTGAGTACATCTGCAGATTTCTGCATGAAATCCAGCTGAAGTGTG
 CCGTCAACAAGATGAGCGTGGACAACCTGGCCACTGTGATCGGGTGAACCTCATCAGGTGGAAGTTGA
 AGACCCAGCTGTGATTATGAGAGGGACTCCTCAGATCAAAGAGTGTGACCATGATGATCAGAGACCAC
 GAAGTCTCTTCCCAAGTCTAAGGATGCACCGATCTCACCCCTGCCAGAAAAACGATGCCAAGAAGG
 CTCCAGTGCCCGAAGCTCTGTGGGTGGGATGCCACGGAAGACCCACCCCTTTCTAGGACAGCAGCTT
 CAGTAACACAGCAAGTAGTCTGATGCCACAGTCCCACTGGACCACTGCCAGTACCAGCATCAGGAA
 GACAGCGGAAAAGCCCCAGGAAAACCCAGGAGACTGGAAGATGCAATCCCGTAAAAGGACTCAAACGC
 TTCCCAACCGGAAGTGTTCCTGACGTCGCATTCCAAGGCACCACCAGCAGTAAACTGGAATCTTTAA
 AAATGAGTTCTGGTCTCCATCTCAGAGGCTAAGGCAGGAGAAGGGCACAGGCGAACTATGTCCCAAGAC
 TTGCGCCACCTTTCCAATGACCAGCGGACTTCTACCTACGATAATGTCCCACCTCACCAGTCCCAAG
 GGAACCCAGCAGGTGCACTCTCTCCCCTGCCAGTACTCAAGAGAGATGCTCTTGTAGCACAGACTC
 TGAATGGAGGCTGGAAGCAAGAACTCTGGCGAGGATGACCTTGATTCTCTGCAGAGGACAGTCCAGAGC
 CTACAGAAGGAAATAGAAACCCAGAAGCAGGTCTATGAGGAACAGATTAACAACTGGAGAGGAAATTT
 ACGATGTCTGGGCTAAGGTGGTGAAGCTCAATGAAGAAGTTCGAGAGGGAGAGGAAGAAATTCGCGGCCCT
 GGAAATCAGCCTTCGAAATGTGGAGCGCTCCCGGAGGACGTTGAGAAGGGAACAGAGTCTTGGAAGAA
 GAAGTCAAGGAGTTTGTGAAGTCGATGGAGAAGCCCAAGACAAAGACGGATCCT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR216611 representing NM_001037727
Red=Cloning site Green=Tags(s)


```
MSLKLPRNWFNLKAEASKIARSRVMTGEQMAAFHPPTTPNPLERPIKVGWLKKQRSIVKNWQQRVFL
RAQQLYYYKDEEDSKPQGCMYLPGSTVKEIATNPEEAGKFVFEVIPASSDQNRIGQDSYVLMASQVEME
EWWKFLRRVAGTPSGAVFGQRLDETVAYEQKFGPHLVPILVEKCAEFILEHGVSEEGIFRLPGQDNLVKQ
LRDAFDAGERPSFDRD TDVHTVASLLKLYLRDLPEPVPVWSQYEGFLLCGQLMNADEAKAQQLVKQLST
LPRDNYNLLSYICRFLHEIQLNCAVNKMSVDNLATVIGVNLIRSKVEDPAVIMRGTPQIQRVMTMMIRDH
EVLFPKSKDAPI SPPAQKNDAKKAPVPRSSVGDATEDPPLSRTDSFSNTASSPDATSPTGPLPSDQHQE
DSGKAPRENPGDWKMQRKRTQTLPNRKCFLTSAFQGTSSKLEIFKNEFWSPSSEAKAGEGHRRTMSQD
LRHLSNDQRTSTYDNVPTSPQSQGNPAGAL SPPASDSKRDALVSTDSEMEAGSKNSGEDDLDSLQRTVQS
LQKEIETQKQVYEEQIKNLEKENYDVWAKVVR LNEELERERKKFAALEISLRNVERSREDVEKRNRLVEE
EVKEFVKSMKPKTKTDP
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001037727

ORF Size: 1944 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001037727.2](#), [NP_001032816.1](#)

RefSeq Size: 2831 bp

RefSeq ORF: 1947 bp

Locus ID: 232201


UniProt ID: [Q8BYW1](#)

Cytogenetics: 6 D1

MW: 73.8 kDa

Gene Summary: GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR216611