

Product datasheet for **MR216611L3V**

Arhgap25 (NM_001037727) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Arhgap25 (NM_001037727) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Arhgap25
Synonyms:	A130039I20Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001037727
ORF Size:	1944 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216611).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001037727.1 , NP_001032816.1
RefSeq Size:	2831 bp
RefSeq ORF:	1947 bp
Locus ID:	232201
UniProt ID:	Q8BYW1
Cytogenetics:	6 D1
Gene Summary:	GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state.[UniProtKB/Swiss-Prot Function]



[View online »](#)