

Product datasheet for **MR216627**

Adgre4 (NM_139138) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Adgre4 (NM_139138) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Adgre4
Synonyms:	D17Ert479e; Egf-tm7; Emr4; Fire; Gpr127
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>MR216627 representing NM_139138
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGTTGATGGGAGCAACTAGAGATATGGGAAGCAGGTGCCTTCTGCATGCCTCAGTTCCTGGAATGCTGC
 TTATCTGGTCAATATTACAAATGATGAATATTTACGTTCTGTCCCCAGTGCAATGAAAATGCCAGCTG
 CTTCAACAGCACCCACTGTGTTTGTAAAGAAGGATTCTGGACGGGCTCTGAGAATAGAAGAATTATTGAG
 CCCCATGAGAAATGTCAAGATATTAATGAGTGTCTACTGAAAGAATTGGTATGCAAGGATGTGCTGACT
 GCAGAAATAAAATTGGGACTTACATATGCAGCTGTGTAGTAAATATCCTTTGTTCAACTGGGTAGCTGG
 CATTATTAATATTGATCACCTGATTGTTATGTGAACAAGAGCAAGAATACAGGATCAAAAACACATACT
 TTGGGAGTACTGAGTGAATTTAAATCAAAGAGGAGGTGCAAAAGGAGCTACCAAGTTACTTCGCAAAAG
 TGAACATCACATCTTGAATGAAAACCTCAGATATACAAAAAAGGATGAAAATCCTTTATTGGATATAGT
 GTATGAACTAAGAGGTGCAAGACGATGACTCTTCTAGAAGCTGGCAACAACACAATGAAGTTGACTGC
 ACTAGTGGTTTCAAAGAGCACAACAGTGGAGGTGAACTGCAGTGGCTTTTCATTGCATATAAGTCTCTTG
 GGAATCTTCTAAATGGTTCCTTTTTTAGTAATGAAGAAGGGTTTCAGGAAGTGACACTGAACTCTCACAT
 CGTTAGTGGAGCCATTGCTCAGAGGTCAAACCTGTCTCTGTAACCTGTACTCCTGACTTTACAAAAT
 ATTCAGCCCATGACTCAAGAGCAGAACATCTCTGTGTCCATTGGGAAGGATCAGAGGAAGGGGGAGCT
 GGTCTACCAAAGGATGCTCTCACGTGTACACCAATAATCCTACACCATTTGCAAGTGTTCCACCTGTC
 CAGCTTTGCTGTGCTCATGGCTCTACCCATGAGGAGGATGGTGTGCTTTCTGCACTCTCTGTGATCACC
 TATGTGGGACTGAGTCTTCTCTCTGTGCCTATTCTGGCGGCCATCACTTTTCTCTGTGCCGACCCA
 TTCAGAATACCAGCAGCACTCCACCTGCAGCTCTCCATCTGCCTTTTCTGGCTGACCTCCTCTTCTCT
 CACAGGCATCAACAGAACTAAGCCTAAGGTGCTGTGCTCCATCATAGCGGGGATGTTGCACTACCTCTAC
 TTGGCTTCCTTCATGTGGATGTTTCTGGAAGGGCTACATCTTTTTCTCACTGTGAGCAATCTCAAAGTGG
 CCAACTACAGCAACTCAGGCAGATTCAAGAAGAGGTTCATGTATCCTGTAGGATATGGGCTTCCTGCTTT
 TATTGTTGCTGTATCTGCAATAGCTGGCCACAAGAATTATGGAACACACAACCACTGCTGGCTCAGCCTT
 CATCGAGGATTCATCTGGAGCTTCTGGGGCCAGCGGCAGCCATTATCTTGATAAACCTGGTGTTCTACT
 TTCTAATAATATGGATTTTGAGAAGCAAACTTTCTCTCTCAATAAAGAAGTTTCTACACTTCAAGACAC
 AAAGGTTATGACATTTAAAGCCATTGTCCAGTTATTTGTGTTGGGATGTTCTTGGGGCATTGGCTTGTTT
 ATTTTCATTGAAGTTGGGAAGACAGTGAGACTGATCGTTGCCTATCTGTTCAACCATCATCAATGTCCTGC
 AGGGTGTTTTGATATTTATGGTACATTGTCTGCTTAATCGCCAGGTGCGGATGGAATATAAGAAGTGGTT
 TCATAGACTGCGGAAGGAAGTTGAAAGTAAAGCACTGAAGTGTCTATTCTACTACTCACACAAAAATG
 GGTCTTTCTCTGAACCTGGAAAATTTCTGCCAACAGGAAACCTCCATGATCCTTCTGACTCCATCCTTC
 CAAGTACTGAAGTAGCAGGTGTATCTAAGCACACCCAGGTCTCACATGGGTGCTGAGGATGTGAACCTC
 AGGTACTCACGCTTACTGGAGCAGAACTATTAGTGAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR216627 representing NM_139138
 Red=Cloning site Green=Tags(s)


MLMGATRDMGSRCLLHASVPGMLLIWSILQMMNISASCPQCENENASCFNSTHCVCKEGFWTGTSENRRRIIE
 PHEKCDINECLLKLCKDVSYCRNKIGTYICSVVYPLFNWVAGIINIDHPCYVYVNSKNTGSKTHT
 LGVLSEFKSKEEVAKGATKLLRKVEHHILNENS DIPKKDENPLLDIVYETKRCKTMTLLEAGNNTMKVDC
 TSGFKEHNSGGETAFAFIAYKSLGNLLNGSFFSNEEGFQEVTLNSHIVSGAIRSEVKPVLSEPVLLTLQN
 IQPIDSRAEHLCVHWEGSEEGSWSTKGC SHVYTNNSTYICKCFHLSSFAVLMALPHEEDGVL SALSVIT
 YVGLSLSLCLFLAAITFLLCRPIQNTSTTLHLQLSICLFLADLLFLTGINRTKPKVLC SIIAGMLHYLY
 LASFMWMFLEGLHLFLT VSNLKVANYSNSGRFKRFMYPGVGLPAFIVAVSAIAGHKNYGTHNHCWLSL
 HRGFIWSFLGPAAAIILINLVFYFLIIWILRSKLSLNKEVSTLQDTKVMTFKAIVQLFVLGCSWGIGLF
 IFIEVGKTVRLIVAYLFTIINVLQGVLI FMVHCLLNQVRMEYKKWFHRLRKEVESESEVSHSTHTTKM
 GLSLNLENFCPTGNLHDPDSILPSTEVAGVYLSTPRSHMGAEDVNSGTHAYWSRTISD

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9013_g12.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_139138


ORF Size: 2067 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_139138.3, NP_631877.2</u>
RefSeq Size:	3491 bp
RefSeq ORF:	2070 bp
Locus ID:	52614
UniProt ID:	<u>Q91ZE5</u>
Cytogenetics:	17 28.98 cM
MW:	77.5 kDa
Gene Summary:	May mediate the cellular interaction between myeloid cells and B-cells.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR216627