

## Product datasheet for **MR216650L3V**

### **Rusc1 (NM\_001083807) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Rusc1 (NM_001083807) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rusc1
Synonyms:	2210403N08Rik; AA408288; NESCA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001083807
ORF Size:	2679 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216650).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001083807.1</a> , <a href="#">NP_001077276.1</a>
RefSeq Size:	3388 bp
RefSeq ORF:	2682 bp
Locus ID:	72296
UniProt ID:	<a href="#">Q8BG26</a>
Cytogenetics:	3 F1



[View online »](#)

**Gene Summary:**

Putative signaling adapter which may play a role in neuronal differentiation. May be involved in regulation of NGF-dependent neurite outgrowth. Proposed to play a role in neuronal vesicular trafficking, specifically involving pre-synaptic membrane proteins. Seems to be involved in signaling pathways that are regulated by the prolonged activation of MAPK. Can regulate the polyubiquitination of IKBKG and thus may be involved in regulation of the NF-kappa-B pathway.[UniProtKB/Swiss-Prot Function]