

Product datasheet for **MR216687L4V**

Nr2c2ap (NM_001025587) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nr2c2ap (NM_001025587) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nr2c2ap
Synonyms:	2310073E15Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001025587
ORF Size:	348 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216687).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001025587.1 , NP_001020758.1
RefSeq Size:	1182 bp
RefSeq ORF:	351 bp
Locus ID:	75692
UniProt ID:	Q3TV70
Cytogenetics:	8 B3.3
Gene Summary:	May act as a repressor of NR2C2-mediated transactivation by suppressing the binding between NR2C2/TR4 and the TR4-response element in target genes.[UniProtKB/Swiss-Prot Function]



[View online »](#)