

Product datasheet for MR216717

Gpr158 (NM_001004761) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Gpr158 (NM_001004761) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Gpr158
Synonyms: 5330427M13Rik; mKIAA1136
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR216717 representing NM_001004761
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGGGAGCCATGGCTTACTCTTTACTCTTCTGCCTCCTGCTCGCTCATCTGGGATTGGGAGAGGTTGGTG
 CCAGCCTCGACCCTCCGGGACGGCCAGATTCCCCTCGAGAGAGGACTCCGAGGGGAAAGCAGCATGGTCA
 GCAGCTGCCGCGAGCCTCTGCCCCAGACCCTCTATTCCCTGGAGCCGCTCCACTGACGGCACCATCTTG
 GCACAAAACTAGCCGAGGAGGTGCCCGTGGACGTGGCCTTTACCTCTACACTGGGACTTCCACCAGC
 TGAAGAGAGCCAACTGCTCAGGCCGCTACGAGTTGGCTGGCCTGCCGGGAAAGTACCGTCCCTAGCCAG
 CTCCCATCCTTCCCTTGCACGGGGCGCTGGACACTCTGACACACGCCACCAACTTCCCTCAACATGATGTTG
 CAGAGCAATAAGTCTCGGGAACAGACCCTGCAGGATGACCTGCAGTGGTACCAGGCTCTGGTGCGAAGTC
 TCCTGGAGGGTGAGCCAGCATCTCCCGGGCGGCCATCACCTTCAGCACCAGAGTCGCTGTCCACGCCCGC
 CCCGCAGGTCTTCTCCAGGCCACCCGCGAGGAGAGTCGATCCTGCTGCAGGACCTGTCTCTCCGCT
 CACCACCTGGCTAACGCCACTCTAGAGACCGAGTGGTTCATGGGCTTCGGCGAAAGTGGAGGCCCACT
 TGCACCGCCGTGGCTCCAATCAGGGGCCACGGGGCCTAGGCCACAGCTGGCGGCGAGGGATGGGCTTGG
 AGGGGACAGGAGCCACGTCAAGTGGTCCCACCTTATCTGGAGTCCGAGAATGGGAGTTACAAGCCTGGG
 TGGCTGGTCACGCTCTCTGCTGCTTCTACGGGCTTCAACCTAACCTGGTCCCAGGAAATCAGGGGTGTA
 TGAAGTTGATATAAGCCTTCAGAAAGTGGACATTGACCAATGTTTCGAGTATGGCTGGTTTTTCAGGAAC
 TCACAAATGCCACCTTAACAACCTCAGAGTGTATGCCAATCAAAGCCTGGGATTCTGCTTGGAGCCTAT
 CAGTGCATTTGCAAAGCAGGGTTCTATCATCCCGTGTCTTCTCAGTGAACAACCTCCAGAGAAGGGGTC
 CAGATCACCATTTTTAGGAAGCACAAAGGATGTGTCAGAAGAAACCCACGTCTGCCTGCCTGCAGAGA
 AGGCTGCCCTTCTGTGCTGATGATCGTCCATGCTTTGTCCAGGAGGATAAATACCTGCGGCTCGCTATT
 ATCTCTTTCCAAGCCCTGTGTATGCTGCTGGACTTTGTCCAGCATGCTGGTGGTCTACCATTTTCGAAAG
 CAAAGAGCATCAGGCATCTGGGCTCATCTGTTGAAACAATATTGTTGGTCTCTGCTCCTATACTT
 TCCGGTCTTATTCTGTACTTTGAGCCAAGTACATTTCTGTGATTCTCCTAAGGTGGGCTCGTCTCTT



GGTTTTGCTACTGTTTATGGAAGTGTGACCCTCAAGCTTCACAGGGTTTTGAAGGTGTTCTTTACGAA
CAGCACAAGAATCCCATACATGACTGGTGGTCGGGTCATGAGGATGCTGGCAGTAATAGTCTTGGTTGT
GTTTTGGTTTCTGTTGGCTGGACTTCTTCTATGTGCCAGAAATTTGGAGCGGGATATTTTCTTGTGGC
CAAGGGCAAACATCTGATCACCTCACCTCAATATGTGCCTCATTGACCCTGGGATTACATGACAGCAG
TTGCTGAATTTTTATTCTCTTGTGGGGCATTATCTCTGCTATGCAGTGCAGTGTCCCTTACAGATT
TCATGAGCCACGCTACATGGCTGTCGAGTTCACAATGAGCTCATTATCACGGCTATATTCACATAAATT
AGATTCGTGCTCGCCTCGAGACTTCAGCCCGACTGGATGCTGATGCTGTACTTTGCACATGCTCACTTGA
CTGTGACAGTCACCATTGGGCTGCTTCTAATTCCAAAGTTTTACATTCAAGCAATAATCCCGAGATGA
CATTGCTACGGAAGCATATGAAGATGAGCTGGACATGGGTCGCTCTGGATCCTACCTGAACAGCAGTATC
AATTCAGCCTGGAGTGAACACAGCCTAGATCCAGAAGACATTCGGGATGAACTGAAAAAAGTCTATGCTC
AGTTAGAAATACAAAACGTAAGAAAGTATGATCACAACAACCCACCTTCAGAAAAACGGTGTCCAA
GAAGGGCTTAGGCCGGTCCATCATGAGGCGCATCACTGAGATCCCAGAGACAGTCAAGCAGCAGTGTCC
AAAGAGGACAAGGAAGGCACAGACCACAGTGCAGCCAAGGGGACAGGCCCTGTCCGGAAGAACCCTACTG
AGTCTCAGGGAATACAGGGAGACCCAAGGAAGAGTCCCTTAAAAATCGGGTCTTCTCTCAAGAAATC
CCACAGCACATATGACCATGTGAGAGACCAAAGTGTAGAAATCCAGCAGCTTGCCCATGGAGACCAGGAG
GAGGAGGCCACCGAAAATCAACCTTGGAGTCGCTGTCAAGTAAAAAAGTGCACAAAAAGCTCAAAGAAG
ATAGTGAGGCCGAGTCCACTGAGTCAGTCCGTTGGTCTGCAAGTCAGCCAGTGTCCACAACCTCAGCTC
AGAGAAGAAGCCTGGGCACCCACGTACATCTATGTTACAGAAATCTCTCAGTGTCAATGCAAGTGCCAAG
GAGAAGACACTTGGCCTGGCTGGCAAAACCCAGACTTTAGTTATGGAAGACCGAGCTAAGTCCCAGAAAC
CAAAAGATAGAGAGACCATCAGGAAGTACTCAAATTCAGATAATGTAGAGACTATCCCAAAGTCCGGTCA
CATGGAGGAGCCAAGAAAGCCCCAGAAATCTGGGATTATGAAACAGCAGAGGGTCAAGCTCCCACTGCC
AATCCTGACGTGAGTCCGGTATCACCCAGATAAAGGACAATTTGATATCGGAGAAGTATGCCCTGGG
AGGTTTATGACCTGACCCCTGGTCTATGCCTTCAGAGCCAAAAGCTCAAAAACACGTATCAATCGCAGC
TTCTGAAGTGGAACAAAACCCCTGCTTCGTTTTTGAAGGAGAAGTCTACCACAAGTCTAAGGCAACTGAA
GGTCTCTACCAAGCCAATCACAAGAGCATAGACAAGACAGAGGTGTGCCCTGGGAAATCCACAGTCAGT
CCCTTTTAGAAGATGAGAATCGTCTGATTTCTAAGACTCCAGTCTCCAGGGAGAGCCAGAGAGGAGAA
TGGAAAGCCAGCTTACACAACCAATATGTGTGCTGGGCAATATGAAGAAGTGCACCACAAAAGTCGTAGCA
CCAAAAGTAGAGAATGAAAATCTCAACCAATGGGAGACCAGGAGAAACAGACATCTTCTCTGTGGATA
TCATTCTGGCTCCTGTAAGTCAAGTAATAACTCCCATCAGCCCTAACATCACGAGCCGAAGTGTGTCC
CTGGGAGTTCGAGCCCTAGAACAACCAATGCTGAAAGAAGCGTAACTTACCTGCCTCCTCTGCTTTA
AGTGCAAATAAGATACCAGGGCCTCAGAAA

ACGCGTACGCGGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR216717 representing NM_001004761
 Red=Cloning site Green=Tags(s)

```

MGAMAYSLLFCLLLAHLGLGEVGASLDPPGRPDSPRERTPRGKQHGQQLPRASAPDPSIPWSRSTDGTIL
AQKLAEEVPVDVASYLYTGDFHQLKRANCSGRYELAGLPGKSPSLASSHPSLHGALDTLTHATNFLNMML
QSNKSREQTVQDDLQWYQALVRSLLLEGEPISIRAAITFSTESLSTPAPQVFLQATREESRILLQDLSSSA
HHLANATLETWFHGLRRKWRPHLHRRGSNQGPRGLGHSWRRRDGLGGDRSHVKWSPPYLECENGSYKPG
WLVTLSSAAFYGLQPNLVPEFRGMKVDISLQKVIDIDQCSSDGFSGTHKCHLNNSECMPIKGLGFVLGAY
QCICKAGFYHPRVFSVNNFQRRGPDHFFSGSTKDVSEETHVCLPCREGCPFCADDRPCFVQEDKYLR LAI
ISFQALCMLLDFVSM LVVYHFRKAKSIRASGLILLETILFGSLLL YFPVVILYFEPSTFR CILLRWARLL
GFATVYGTVT LKLRVLKVFLSRTAQRIPYMTGGRVMRMLAVIVLVVFWFLVGWTS SMCQNLERDILLVG
QGQTS D H L T F N M C L I D R W D Y M T A V A E F L L W G I Y L C Y A V R T V P S A F H E P R Y M A V A V H N E L I I T A I F H T I
R F V L A S R L Q P D W M L M L Y F A H A H L T V T V T I G L L L I P K F S H S N N P R D D I A T E A Y E D E L D M G R S G S Y L N S S I
N S A W S E H S L D P E D I R D E L K K L Y A Q L E I Y K R K K M I T N N P H L Q K K R C S K K G L G R S I M R R I T E I P E T V S R Q C S
K E D K E G T D H S A A K G T G L V R K N P T E S S G N T G R P K E E S L K N R V F S L K K S H S T Y D H V R D Q T V E S S L P M E S Q E
E E A T E N S T L E S L S S K K L T Q K L K E D S E A E S T E S V P L V C K S A S A H N L S S E K K P G H P R T S M L Q K S L S V I A S A K
E K T L G L A G K T Q T L V M E D R A K S Q K P K D R E T I R K Y S N S D N V E T I P N S G H M E E P R K P Q K S G I M K Q R V S L P T A
N P D V S S G I T Q I K N F D I G E V C P W E V Y D L T P G P M P S E P K A Q K H V S I A A S E V E Q N P A S F L K E K S Y H K S K A T E
G L Y Q A N H K S I D K T E V C P W E I H S Q S L L E D E N R L I S K T P V L P G R A R E E N G S Q L Y T T N M C A G Q Y E E L P H K V V A
P K V E N E N L N Q M G D Q E K Q T S S S V D I I P G S C N S S N N S H Q P L T S R A E V C P W E F E P L E Q P N A E R S V T L P A S S A L
S A N K I P G P Q K
  
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9001_g09.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:




ACCN: NM_001004761

ORF Size: 3600 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001004761.1 , NP_001004761.1
RefSeq Size:	7143 bp
RefSeq ORF:	3603 bp
Locus ID:	241263
UniProt ID:	Q8C419
Cytogenetics:	2 A3
MW:	134.9 kDa
Gene Summary:	Orphan receptor.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR216717