

Product datasheet for **MR216730L4V**

Taf7l (NM_028958) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Taf7l (NM_028958) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Taf7l
Synonyms:	4933438I11Rik; Taf2q
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_028958
ORF Size:	1413 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216730).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_028958.4 , NP_083234.2
RefSeq Size:	1935 bp
RefSeq ORF:	1416 bp
Locus ID:	74469
UniProt ID:	Q9D3R9
Cytogenetics:	X E3



[View online »](#)

Gene Summary:

Probably functions as a spermatogenesis-specific component of the DNA-binding general transcription factor complex TFIID, a multimeric protein complex that plays a central role in mediating promoter responses to various activators and repressors. May play a role in spermatogenesis.[UniProtKB/Swiss-Prot Function]