

## Product datasheet for **MR216775L3V**

### Ap3m2 (NM\_001122820) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ap3m2 (NM_001122820) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ap3m2
Synonyms:	5830445E16Rik; AP-3B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001122820
ORF Size:	1254 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216775).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001122820.1</a> , <a href="#">NP_001116292.1</a>
RefSeq Size:	3467 bp
RefSeq ORF:	1257 bp
Locus ID:	64933
UniProt ID:	<a href="#">Q8R2R9</a>
Cytogenetics:	8 A2



[View online »](#)

**Gene Summary:**

Component of the adaptor complexes which link clathrin to receptors in coated vesicles. Clathrin-associated protein complexes are believed to interact with the cytoplasmic tails of membrane proteins, leading to their selection and concentration. Ap47 is a subunit of the plasma membrane adaptor (By similarity). In concert with the BLOC-1 complex, AP-3 is required to target cargos into vesicles assembled at cell bodies for delivery into neurites and nerve terminals.[UniProtKB/Swiss-Prot Function]