

Product datasheet for **MR216810A1V**

Mouse Mrfp1 (NM_026242) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Mrfp1 (NM_026242) AAV Particle
Tag: Myc-DDK
Symbol: Mrfp1
Synonyms: 9130413I22Rik; AL022827; Pam14
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR216810 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGCGGCCCTGGACGCGGTGGAGCTGGCGGAGCCCGAGGAGGTGGAGGTGCTGGAGCCCGAGGAGGACT
 TCGAGCAGTTTCTGCTGCCCGTCATCCACGAGATGCGCGAGGACATCGCGTCGCTGACGCGCGAGCGCGG
 GCGCGCGCCGGCGCGCAACCGGGGAAGCTGTGGGAGATGGACAATATGCTGATCCAGATCAAGACGCAG
 GTCGAGGCCTCCGAGGAGAGCGCCCTCAACCACCTGCAGGGTGCCGGCGGCGCCGAGCCCGCGGCCCC
 GGGCGGAGAAGGCCGACGAGAAGGCGCAGGAGATGGCGAAGATGGCCGAGATGCTGGTGCAGCTCGTGCG
 GCGGATAGAGAAGAGCGAGTCTTCG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR216810 protein sequence
 Red=Cloning site Green=Tags(s)

MRPLDAVELAEPEEEVEVLEPEEDFEQFLLPVIHEMREDIASLTRERGRAPARNRGKLWEMDNMLIQIKTQ
 VEASEESALNHLQGAGGAEPGRPRAEKADEKAQEMAKMAEMLVQLVRRIEKSESS

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_026242


[View online »](#)

ORF Size:	375 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_026242.3, NP_080518.1</u>
RefSeq Size:	1630 bp
RefSeq ORF:	378 bp
Locus ID:	67568
UniProt ID:	<u>Q9CQL7</u>
Cytogenetics:	5 B3
MW:	14.2 kDa