

## Product datasheet for **MR216836L3V**

### Mgat1 (NM\_001110148) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Mgat1 (NM_001110148) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mgat1
Synonyms:	Mgat-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001110148
ORF Size:	1341 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216836).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001110148.1</a> , <a href="#">NP_001103618.1</a>
RefSeq Size:	3165 bp
RefSeq ORF:	1344 bp
Locus ID:	17308
UniProt ID:	<a href="#">P27808</a>
Cytogenetics:	11 29.27 cM
Gene Summary:	Initiates complex N-linked carbohydrate formation. Essential for the conversion of high-mannose to hybrid and complex N-glycans.[UniProtKB/Swiss-Prot Function]



[View online »](#)