

Product datasheet for MR216842L4

Prkd3 (NM_001171005) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Prkd3 (NM_001171005) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Prkd3
Synonyms:	4930557O20Rik; 5730497N19Rik; Pkcnu; PKD3; Prkcn
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216842).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	NM_001171005
ORF Size:	2385 bp



[View online »](#)

