

## Product datasheet for **MR216845L3V**

### Hagh (NM\_024284) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Hagh (NM_024284) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hagh
Synonyms:	BC019817; Glo-2; Glo2; Rsp29
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_024284
ORF Size:	927 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216845).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_024284.2</a>
RefSeq Size:	1197 bp
RefSeq ORF:	930 bp
Locus ID:	14651
UniProt ID:	<a href="#">Q99KB8</a>
Cytogenetics:	17 12.53 cM
Gene Summary:	Thiolesterase that catalyzes the hydrolysis of S-D-lactoyl-glutathione to form glutathione and D-lactic acid.[UniProtKB/Swiss-Prot Function]



[View online »](#)