

Product datasheet for **MR216861L3V**

Klhl40 (NM_028202) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Klhl40 (NM_028202) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Klhl40
Synonyms:	2310024D23Rik; Kbtbd5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_028202
ORF Size:	1863 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216861).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_028202.3 , NP_082478.1
RefSeq Size:	2455 bp
RefSeq ORF:	1866 bp
Locus ID:	72330
UniProt ID:	Q9D783
Cytogenetics:	9 F4


[View online »](#)

Gene Summary:

Substrate-specific adapter of a BCR (BTB-CUL3-RBX1) E3 ubiquitin ligase complex that acts as a key regulator of skeletal muscle development (PubMed:25940086). The BCR(KLHL40) complex acts by mediating ubiquitination and degradation of TFDP1, thereby regulating the activity of the E2F:DP transcription factor complex (PubMed:25940086). Promotes stabilization of LMOD3 by acting as a negative regulator of LMOD3 ubiquitination; the molecular process by which it negatively regulates ubiquitination of LMOD3 is however unclear (PubMed:24960163).[UniProtKB/Swiss-Prot Function]