

Product datasheet for **MR216923L4V**

PPP1R1C (NM_172420) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PPP1R1C (NM_172420) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	PPP1R1C
Synonyms:	4930565M23Rik; 4930579P15Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_172420
ORF Size:	327 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR216923).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_172420.3 , NP_766008.2
RefSeq Size:	3071 bp
RefSeq ORF:	327 bp
Locus ID:	75276
UniProt ID:	Q8BKK4
Cytogenetics:	2 C3
Gene Summary:	May increase cell susceptibility to TNF-induced apoptosis.[UniProtKB/Swiss-Prot Function]



[View online »](#)